



E100 Learning Wide Field Microscope Halogen Illumination Microscopes A12.0705

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.0705
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Trinocular
- Illumination: Halogen Lamp
- Product Name: Compound Optical Microscope
- Highlight: **binocular compound microscope ,
optical compound microscope**



for more products please visit us on cnoec.com

Product Description

E100 Learning Wide Field Microscope Halogen Illumination Microscopes A12.0705


A12.0705 Biological Microscope	
Optical System	CFI Optical System
Head	E2-TB Binocular Viewing Head, Inclined at 30°, Interpupillary 47-75mm
Eye-piece	E1-CFI W10x/18
Objective	CFI BE Plan Achromatic Objective 4x ,10x , 40x , 100x
Nosepiece	Frontward Quadruple Nosepiece
Focusing	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke 37.7mm per Rotation, Fine Stroke 0.2mm per Rotation, Moving Range 22mm
Stage	Double Layers Mechanical Stage 155*134mm
Condenser	Abbe Condenser NA1.25 can hold 33mm filter
Illumination	6V/20W Halogen Lamp
Optional Accessories	
	Item No.
Head	E2-TF Trinocular Viewing Head, Inclined at 30°, 47-75mm
Eye-piece	CFI E 15x/12
Objective	ICFI BE Plan Achromatic Objective 20X , 60X
Phase Contrast	Phase Contrast Kit
Pointer	
Condenser	Dark Field Condensor
Illumination	Reflected Mirror
Package	Wooden box or aluminium alloy Box



A12.0705 Laboratory Biological Microscope



CFI BE Plan Achromatic Objective
4x, 10x, 40x, 100x



Easy Replacment of 6V/20W
Halogen Lamp



Double Layers Mechanical Stage
155*134mm



A12.0705 Laboratory Biological Microscope



A12.0705 Laboratory Biological Microscope



CFI BE Plan Achromatic Objective 4x,10x,40x,100x



Easy Replacment of 6V/20W Halogen Lamp



Double Layers Mechanical Stage 155*134mm



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China