



Lab Achromatic LED Biological Microscope Monocular Microscopes A11.1523

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1513
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student
- Certification: CE
- Drawtube: Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **binocular biological microscope,
biological compound microscope**



More Images



Product Description

Lab Achromatic LED Biological Microscope Monocular Microscopes A11.1523

A11.1523 Student Microscope, Dual LED Light Source, New developed model in 2014, an ideal tool or gift scope for children and student, with full accessories kit all in one..

- Full Metal Body With High Quality Optical Glass Makes It A Real Microscope
- Dual LED Light Source, Top and Bottom, Brightness Adjustable
- Powered by 3A Battery, Cordless For Easy Moving And Outdoor Use
- Order Shipped By Air Will Supply WITHOUT Battery

A11.1523 Student Microscope, Dual LED Light Source	
Head	Monocular Head with Metal Eyepiece Tube Dia. 23.2mm, 360° Rotatable
Eyepiece	WF10x
Objective	Achromatic 4x,10x,40x Spring
Nosepiece	Triple
Working Stage	Metal Plain Stage, With Clips, Size 90*90mm
Condenser	Disc Diaphragm, 0.65 Condenser
Illumination	Dual LED Light Source, Top & Bottom, Brightness Adjustable
Power	3A Batteries
Accessories Kits	--Lens Paper, --Dissecting Needle, --Cotton Stick, --Plastic Petri Dish, --Plastic Test Tube, --Plastic Forceps, --Prepared Slide 5 pcs, --Blank Slide 5 pcs, --Cover Glass 5 pcs, --Slide Sticker 5 pcs, --Red Stain (food color), -- Blue Stain (food color) --Manual Microtome



A11.1523 Student Microscope. Dual LED Light Source



A11.1523 Student Microscope, Dual LED Light Source, With Handle For Easy Carry



