



200X 1.3M USB Handheld Digital Microscope Digital Camera Microscopes A34.5003

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.5003
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE, Rohs
- Magnification: 200x
- Theory: Handheld Digital Microscope
- Highlight: **Hand held Digital Microscope, hand held microscope**



More Images



Product Description

200X 1.3M USB Handheld Digital Microscope Digital Camera Microscopes A34.5003

- A34.5003 USB Digital Microscope, 200X, 1.3M
- Magnification 20X~200X
- Resolution 1.3M, Upto 2.0M By Software
- 8 LED, Adjustable By Control Wheel
- Manual Focus From 10mm To 500mm

Our Product

A34.5003 USB Digital Microscope	
Image Sensor	1.3 Mega Pixels, (Interpolated to 2M)
Still Capture Resolution	1600*1200,1280*1024,1280*960,1024*768, 800*600,640*480,352*288,320*240,160*120
Video Capture Resolution	1600*1200,1280*1024,1280*960,1024*768, 800*600,640*480,352*288,320*240,160*120
Focus Range	Manual Focus From 10mm To 500mm
Frame Rate	Max 30fps Under 600 Lus Brightness
Magnification Ratio	20x-200x
Video Format	AVI
Photo Format	JPEG / BMP
Light Source	8 LED, Adjustable By Control Wheel
PC Interface	USB2.0
Power Source	5V DC From USB Port
Operation System:	Windows 2000/XP/Vista/Win7
	Mac OS X 10.5 or Above
OSD Language	English, German, Spanish, Korean,French, Russian
Bundle Software	MicroCapture With Measurement
Product Size	110mm (L) x 33mm (R)
Product Weight	160g
Packing List	Microscope Metal Stand Software CD User's Manual
Package	20Pcs/Carton, 10Kgs/Carton, 46(L)*35(W)*30.0(H)cm
Certificate	With CE, FCC Approval, RoHS Compliance







www.cnobec.com



www.cnobec.com





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China