



USB Handheld Digital Microscope 2.0M

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.5004
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE|Rohs
- Magnification: 200x
- Theory: Handheld Digital Microscope
- Highlight: **Portable Digital Microscope, hand held microscope**



More Images



Product Description

USB Handheld Digital Microscope 200X 2.0M , Digital Camera Microscopes A34.5004

- A34.5004 USB Digital Microscope, 200X, 2.0M
- Magnification 20X~200X
 - Resolution 2.0M, Upto 5.0M By Software
 - 8 LED, Adjustable By Control Wheel
 - Manual Focus From 10mm To 500mm

Our Product

A34.5004 USB Digital Microscope	
Image Sensor	2.0 Mega Pixels, (Interpolated to 5.0M)
Still Capture Resolution	2560*2048 (5M), 2000*1600,1600*1200,1280*1024,1280*960, 1024*768, 800*600,640*480,352*288,320*240
Video Capture Resolution	2560*2048 (5M), 2000*1600,1600*1200,1280*1024,1280*960, 1024*768, 800*600,640*480,352*288,320*240
Focus Range	Manual Focus From 10mm To 500mm
Frame Rate	Max 30fps Under 600 Lus Brightness
Magnification Ratio	20x~200x
Video Format	AVI
Photo Format	JPEG / BMP
Light Source	8 LED, Adjustable By Control Wheel
PC Interface	USB2.0
Power Source	5V DC From USB Port
Operation System:	Windows 2000/XP/Vista/Win7
	Mac OS X 10.5 or Above
OSD Language	English, German, Spanish, Korean,French, Russian
Bundle Software	MicroCapture With Measurement
Product Size	110mm (L) x 33mm (R)
Product Weight	160g
Packing List	Microscope Metal Stand Software CD User's Manual
Package	20Pcs/Carton, 10Kgs/Carton, 46(L)*35(W)*30.0(H)cm
Certificate	With CE, FCC Approval, RoHS Compliance



A34.5004 USB Digital Microscope, 200X,2.0M





A34.5004 USB Digital Microscope, 200X, 2.0M, Zoom Knob



A34.5004 USB Digital Microscope, 200X, 2.0M, LED Light Control Knob



A34.5004 USB Digital Microscope, 200X, 2.0M. 8 LED Light Source





A34.5004 USB Digital Microscope, Packing Box



A34.5004 USB Digital Microscope, 200X,2.0M



A34.5004 USB Digital Microscope, 200X, 2.0M



A34.5004 USB Digital Microscope, 200X, 2.0M, Zoom Knob



A34.5004 USB Digital Microscope, 200X, 2.0M, LED Light Control Knob



A34.5004 USB Digital Microscope, 200X, 2.0M



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China