



Binocular Jewelry Microscope With Zoom Ratio 1:7 A24.0402 For Research

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A24.0402
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Research
- Certification: CE
- Drawtube: Binocular
- Theory: Jewelry Microscope
- Highlight: **laboratory binocular microscope, gemological microscope**

for more products please visit us on cnoec.com

Product Description

Binocular Jewelry Microscope With Zoom Ratio 1:7 A24.0402 For Research

A 24.0402 Specification of Gem Microscope	
Head	Binocular head, interpupillary distance: 52-75mm, diopter adjustment: ± 5 diopter
Body	Optic head and stand tilt to a more natural 40° viewing angle
Eyepiece	SWF10x/22mm
Zoom Ratio	1:7
Objective Zoom Range	6.5x-45x
Focusing	Focus range: 135mm
Working Distance	103mm
Illumination	Voltage input 90-230V
	Adjustable iris diaphragm to control width of 2-42mm light beam
	Built-in optical-fiber illumination system
	Tungsten-halogen darkfield/brightfield illumination



Our Product Introduction

for more products please visit us on cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China