



400X 360°Rotatable Inverted Optical Microscope A14.0701 Binocular Head

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A14.0701
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Max Magnification: 400X
- Product Name: Inverted Optical Microscope
- Theory: Inverted Microscope
- Highlight: **inverted fluorescence microscope, inverted light microscope**

for more products please visit us on cnoec.com

Product Description

400X 360°Rotatable Inverted Optical Microscope A14.0701 Binocular Head

A14.0701 Inverted Microscope		
Optical System	Infinite Optical System	
Head	Siedentopf binocular viewing head, inclined at 45°, 360° Rotatable, Interpupillary 48-75mm	
Eyepiece	Φ30 Wide Field Eyepiece WF10x/22, Eyepiece Tube Diameter 30mm	
Objective	Infinite Plan Achromatic Objective 4x/0.10 WD=22	
	Infinite LWD Plan Achromatic Phase Objective 10x/0.25 WD=6	
	Infinite LWD Plan Achromatic Phase Objective 20x/0.4 WD=3.1	
	Infinite LWD Plan Achromatic Objective 40x/0.55 WD=2.2	
Focusing	Coaxial coarse and fine adjustment, fine division 2μm, Moving Range up 4.5mm down 4.5mm	
Nosepiece	Backward Quintuple Nosepiece	
Condenser	ELWD Condenser NA0.3, LWD72mm (without Condenser 150mm)	
Centering Telescope	Centering Telescope (φ30mm)	
Phase Annulus	10x, 20x, 40x Phase Annulus Plate (Adjustable)	
Stage	Plain Stage 230x170	
	Glasse Insert	
	Auxiliary Stage	
	Terasaki Holder	
	Petri Dish Holder φ38	
	Petri Dish Holder φ54	
	Attachable Mechanical Stage, X, Y Coaxial Control, Moving Rang 80x120mm	
Illumination	6V/30W Halogen Lamp, Brightness Adjustable	
Filter	Blue, Green and ground glass, Diameter 45mm	
	Optional Accessories	Item No.
Eyepiece	Wide Field Eyepiece WF15x/16	
	Wide Field Eyepiece WF20x/12	
Objective	Lamp House Adjustment Objective	
Illumination	Epi-Fluorescent Attachment	
Photo attachment		
C-Mount	0.4x	






Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China