



## Inverted Fluorescence Binocular Compound Microscopes 40X - 400X A16.0701

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A16.0701
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Highlight: **inverted fluorescence microscope, epifluorescence microscope**

## Product Description

Inverted Fluorescence Binocular Compound Microscopes 40X - 400X A16.0701

A16.0701 Inverted Fluorescent Microscope	
Optical System	Infinite Optical System
Head	Siedentopf binocular viewing head, inclined at 45°, 360° Rotatable, Interpupillary 48-75mm
Eyepiece	Φ30 Wide Field Eyepiece WF10×/22, Eyepiece Tube Diameter 30mm
Objective	Infinite Plan Achromatic Objective 4×/0.10 WD=22
	Infinite LWD Plan Achromatic Phase Objective 10×/0.25 WD=6
	Infinite LWD Plan Achromatic Phase Objective 20×/0.4 WD=3.1
	Infinite LWD Plan Achromatic Objective 40×/0.55 WD=2.2
Focusing	Coaxial coarse and fine adjustment, fine division 2μm, Moving Range up 4.5mm down 4.5mm
Nosepiece	Backward Quintuple Nosepiece
Condenser	ELWD Condenser NA0.3, LWD72mm (without Condenser 150mm)
Centering Telescope	Centering Telescope (φ30mm)
Phase Annulus	10×, 20×, 40× Phase Annulus Plate (Adjustable)
Stage	Plain Stage 230×170
	Glasse Insert
	Auxiliary Stage
	Terasaki Holder
	Petri Dish Holder φ38
	Petri Dish Holder φ54
	Attachable Mechanical Stage, X, Y Coaxial Control, Moving Rang 80×120mm
Illumination	6V/30W Halogen Lamp, Brightness Adjustable
Fluorescent Attachment	Epi-Fluorescent Illumination
	B, G
	Blue, Green and ground glass, Diameter 45mm



Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China