



## 0.4MP 480 x 480 Digital Eyepiece Camera A59.1008 With 7Inch LCD Screen

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.1008
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Microscope
- Certification: CE
- Product Name: Microscope Accessories
- Highlight: microscope accessory, parts of the microscope



### More Images



## Product Description

### 0.4MP 480 x 480 Digital Eyepiece Camera A59.1008 With 7Inch LCD Screen

**A59.1008 7" LCD Digital Eyepiece Camera**, Designed for all kinds microscope, telescope use, sharing micro image with your friends and colleagues easily.

--7" LCD Screen


<b>A59.1008 7" LED Digital Camera</b>	
Digital Demonstrator	7" LED screen with resolution of 480*480
Adapter	Dia.23.2mm adapter
Resolution	0.4MP (720*576)
<b>Optional Accessories</b>	
Sytem	Windows 2000, Windows XP, Vista, WIN7(32bit)
Software	Scopeimage plus
TV Cable	RCA Cable
Video Capture Card	Easycap DC60



A59.1008 LED Digital Camera



It can be combined with monocular/ binocular/ trinocular microscope







**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China