



Achromatic Biological Compound Microscope Monocular A11.1134

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1134
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student
- Certification: CE
- Drawtube: Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
biological compound microscope**



Product Description

Achromatic Biological Compound Microscope Monocular A11.1134

A11.1134 Student Microscope, New developed model in 2014, an ideal tool or gift scope for children and student, with full accessories kit all in one..

- 40x~400x Magnification
- Dual Up & Bottom LED Light Source
- Cordless Powered by AA 1.5V Battery 3 pcs
- Full Metal Body for Steady Observation
- Compact All in One Design for Easy Moving

A11.1134 Economical Cordless Student Microscope,40x-400x	
Head	Monocular Head, 360° Rotatable
Eyepiece	WF10x
Objective	Achromatic 4x,10x,40x Spring
Nosepiece	Triple Revolving Nosepiece
Stage	Plan Stage With Clips
Condenser	Disc Diaphragm
Focusing	Coarse Focusing Knob
Illumination	Dual Light Top & Bottom LED Light
Power	By AA 1.5V Battery 3 Pcs.
Remark	Order Shipped By Air Will Not Include Batteries



A11.1134 Economical Cordless Student Microscope,40x-400x





A11.1134 Economical Cordless Student Microscope, 40x-400x



A11.1134 Economical Cordless Student Microscope, 40x-400x

www.cnoec.com



A11.1134 Economical Cordless Student Microscope, 40x-400x



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China