



## LED Light Source Biological Microscope Monocular Microscopes A11.1135

### Our Product Introduction

#### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1135
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



#### Product Specification

- Applications: Student
- Certification: CE
- Drawtube: Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,  
biological compound microscope**



## Product Description

### Lab Achromatic Biological Microscope Monocular Microscopes A11.1135

**A11.1135 Student Microscope**, New developed model in 2014, an ideal tool or gift scope for children and student, with full accessories kit all in one..

- 40x~400x Magnification
- LED Light Source
- Full Metal Body for Steady Observation
- Coaxial Coarse & Fine Focusing Knob

A11.1135 Student Biological Microscope	
Head	Monocular Head
Eyepiece	WF10x
Objective	Achromatic 4x,10x,40x Spring
Nosepiece	Backward Triple Revolving Nosepiece
Stage	Plan Stage With Clips
Condenser	Abbe Condenser N.A.1.25
Focusing	Coaxial Coarse & Fine Focusing Knob
Illumination	LED Light Source



A11.1135 Biological Microscope





A11.1135 Biological Microscope



A11.1135 Biological Microscope



A11.1135 Biological Microscope



**Opto-Edu (Beijing) Co., Ltd.**



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China