



Laboratory Compound Optical Microscope Halogen Illumination Microscopes A12.1501

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Binocular|Trinocular
- Illumination: Halogen Lamp
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **binocular compound microscope, wide field microscope**



More Images



Product Description

Laboratory Compound Optical Microscope Halogen Illumination Microscopes A12.1501

A12.1501 Laboratory Biological Microscope

- Infinity Optical System
- Infinite Plan Achromatic Objective 4X,10X,40X,100X
- Abbe Condensor NA1.25 with Aperture Diaphragm
- Backward Quadruple Nosepiece
- 12V/30W Halogen Lamp, Brightness Adjustable

A12.1501 Laboratory Biological Microscope		A12.1501-B	A12.1501-T
Optical System	Infinity Optical System		
Head	Seidentopf Binocular Viewing Head, Inclined at 30°, Interpupillary 48-75mm	o	
	Seidentopf Trinocular Viewing Head, Inclined at 30°, Interpupillary 48-75mm		o
Eyepiece	WF10x/20mm		
Nosepiece	Backward Quadruple Nosepiece		
Objective	Infinite Plan Achromatic Objective 4X,10X,40X,100X		
Condensor	Abbe Condensor NA1.25 with Aperture Diaphragm		
Focusing	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Moving Range 22mm		
Working Stage	Double Layer Mechanical Rectangular Stage 175x145mm, Cross Travel 77x50mm		
Light Source	12V/30W Halogen Lamp, Brightness Adjustable		
A12.1501 Laboratory Biological Microscope Optional Accessories			
Objective	Infinite Plan Achromatic Objective 20X		
	Infinite Plan Achromatic Objective 60X		
Adapter	C-Mount 0.4x		



A12.1501-T Trinocular Laboratory Biological Microscope

A12.1501-B Binocular Laboratory Biological Microscope







A12.1501-T Trinocular
Laboratory Biological Microscope

A12.1501-B Binocular
Laboratory Biological Microscope




A12.1501 Laboratory Biological Microscope,
With Infinity Optical Objectives



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China