



## WIFI 365X 5.0M Hand Held Digital Microscope For iPad / PC / Android A34.4903

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.4903
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Certification: CE
- Highlight: Portable Digital Microscope, hand held microscope



for more products please visit us on cnoec.com

### More Images



## Product Description

### WiFi 365X 5.0M Hand Held Digital Microscope For iPad / PC / Android A34.4903


**A34.4903 WiFi Microscope, 365X,5.0M, iPad/Android/Win**, High quality upon very competitive cost.

- 0.7~5x, Continus Zoom Lens
- Objective Lens 1.0x
- Working Distance 105mm
- 5.0M CMOS High Resolution Camera
- WiFi IEEE 802.11b/g/n
- Frame Rate Upto 40 fps
- Support Microsoft Windows, MAC OS, IOS, Andriod or Above
- Enable Your Microscope Upgrade to WiFi Digital Microscope Easily!Auto/Manual Exposure
- Capturing Software, Free Driver Included

<b>A34.4903 WiFi Microscope, 365X,5.0M, iPad/Android/Win</b>	
Microscope System	Monocular Zoom Video Microscope <b>A21.0902</b> --0.7~5x, Continus Zoom Lens --Objective Lens 1.0x --Working Distance 105mm --Large Black Pole Stand, With Focus Arm --Standard C-Mount --No Light Source
Camera Size	Dia. 60mm * Height 70mm
Sensor	5.0M True Color 1/2.5"CMOS
Resolution	2592 x 1944 ( 5.0M)
Frequency Range	2.4000--2.4835GHz
Working Modes	A.USB mode B.WiFi mode
Focus Mode	Real-time Auto-focus
Exposure	Wide Range Auto/manual Exposure
White Balance	Automatic White Balance
Function	Video Recording/ Audio Recording And Ttransmission, Photo Snapshot Function
Turn on WiFi Devices	Instant Monitoring From Smartphone/ Ttablet
Support Browser	Safari / Chrome / IE
Frame Rate	JPG Output: 640x480@40fps,800x600@40fps,1280x960@40fps, 1280x720@40fps, 2592x1944@15fps
Transmit Bandwidth	150M Bps
WiFi	IEEE 802.11b/g/n
Power supply	Separable and rechargeable lithium battery(3.7V 1900mAh); DC 5V USB power supply by PC External DC 5V 1.5A adapter
Operating System	Microsoft® Windows® XP / Vista / 7 / 8 (32 & 64 bit) MAC OS X (Wi-Fi Working Mode) IOS5/IOS6(Wi-Fi Working Mode) Andriod 2.2 or above
Hardware Requirement	CPU: Intel Core2 1.6GHz or Above Memory:2GB or Above USB Port:2.0 DVD-ROM(Required When Installing Software)



[www.cnoec.com](http://www.cnoec.com)



**A34.4903 WIFI Microscope, 365X,5.0M,  
For iPad/Android/Win**



www.cnoec.com

**A34.4903 WIFI Microscope, Consist of  
A59.4903 WIFI Digital Camera + A21.0902 Mono Video Microscope**





**A59.4903 WIFI Digital Camera**



**365X,5.0M, For iPad/Android/Win**



[www.cnoec.com](http://www.cnoec.com)

**A21.0902 Monocular Video Microscope**  
0.7~5x, Continus Zoom Knob, 1.0x Lens,  
Working Distance 105mm



**A59.4903 WIFI Digital Camera, 5.0M, iPad/Android/Win**  
With 5.0M High Resolution Digital Camera, WIFI Transfer Mode And U Cam Software, This WIFI Digital Camera Can Help You To Preview & Capture Image On Your iPhone, iPad or Android Mobile Phone!

You Can Easily View And Share The Excellent World Together With Your Friends Now!



Back    Auto exposure    Zoom out    Zoom in    Auto focus    Snapshot    Video    Folder  
Decrease exposure    Increase exposure

[www.cnoec.com](http://www.cnoec.com)



**A34.4903 WIFI Microscope, 365X, 5.0M, iPad/Android/Win**



**A34.4903 WIFI Microscope, Consist of  
A59.4903 WIFI Digital Camera + A21.0902 Mono Video Microscope**

A59.4903 WIFI Digital Camera, 5.0M,iPad/Android/Win



www.cnoec.com



A34.4903 WIFI Microscope, Consist of  
A59.4903 WIFI Digital Camera +  
A21.0902 Mono Video Microscope

A59.4903 WIFI Digital Camera, 5.0M,iPad/Android/Win

With 5.0M High Resolution Digital Camera, WIFI Transfer Mode And U Cam Software, This WIFI Digital Camera Can Help You To Preview & Capture Image On Your iPhone/iPad or Android Mobile Phone!

You Can Easily View And Share The Excellent World Together With Your Friends Now!



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China