



USB 365X 5.0M Digital Microscope For PC Handheld Digital Microscope A34.4904

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.4904
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Highlight: Hand held Digital Microscope, hand held microscope



More Images



for more products please visit us on cnoec.com

Product Description

USB 365X 5.0M Digital Microscope For iPad / PC / Android A34.4904

A34.4904 USB Digital Miroscope, including monocular zoom microscope and high resolution digital camera, is ideal tool for inspect all kinds objects in industry producing line, laboratory, or school.

- 0.7~5x, Continus Zoom Lens
- Objective Lens 1.0x
- Working Distance 105mm
- 5.0M CMOS High Resolution Camera
- 720P HD Video Up To 30fps
- True Color Restore
- Auto White Balance
- Auto/Manual Exposure
- Capturing Software, Free Driver Included

| A34.4904 USB Digital Miroscope, 365X,5.0M,720p | |
|---|---|
| Microscope System | Monocular Zoom Video Microscope A21.0902 --0.7~5x, Continus Zoom Lens --Objective Lens 1.0x --Working Distance 105mm --Large Black Pole Stand, With Focus Arm --Standard C-Mount --No Light Source |
| Digital Camera | 5.0M CMOS High Resolution Digital Camera |
| Senser | 1/2.5" CMOS, True Color |
| Resolution | 2592*1944 Max. |
| Pixel Size | 2.2x2.2µm |
| Frame Rate | @640X 480: 30fps; @2592X1944: 15fps 720P HD Video Up To 30fps |
| White Balance | Auto |
| Exposure | Manual/Auto |
| Image Format | JPG/BMP/PNG/GIF |
| Video Format | AVI |
| Software | Capturing Software, Free Driver |
| System | Microsoft Windows XP/2003/Vista/Win7/Win8 |
| Lens Mount | Standard C-mount |



Standard C-Mount
With Simple Installation Design, It Can Be Used in Trinocular Photo Tube With C-Mount Of All Kinds Biological Microscopes Or Stereo Microscopes

A59.4904 Digital Camera, 5.0M, USB2.0



5.0M High Resolution CMOS Sensor
Adopt The New CMOS Ture 5.0M Sensor, It Greatly Reduce Image Noise.





A59.4904-M USB Digital Camera, 5.0M,720p
 With Dia. 23.2mm Ring For Biological Microscope Eyepeice Tube
 With Dia. 30 and 30.5mm Ring For Stereo Microscope Eyepeice Tube



0.5X And 1.0X Lens Are Equipped With Rubber Ring To Make The Adapter Ring Very Stable Avoid Detachment Troubles



A59.4904 Microscope Digital Camera
 5.0M CMOS High Resolution Camera
 720P HD Video Up To 30fps
 True Color Restore
 Auto White Balance
 Auto/Manual Exposure
 Capturing Software, Free Driver Included



A59.4904 Digital Camera, 5.0M, USB2.0




Standard C-Mount
With Simple Installation Design, It Can Be Used in Trinocular Photo Tube With C-Mount Of All Kinds Biological Microscopes Or Stereo Microscopes



5.0M High Resolution CMOS Sensor
Adopt The New CMOS Ture 5.0M Sensor, It Greatly Reduce Image Noise.



A59.4904-M USB Digital Camera, 5.0M, 720p
With Dia. 23.2mm Ring For Biological Microscope Eyepeice Tube
With Dia. 30 and 30.5mm Ring For Stereo Microscope Eyepeice Tube



0.5X And 1.0X Lens Are Equipped With Rubber Ring To Make The Adapter Ring Very Stable Avoid Detachment Troubles

A59.4904 Microscope Digital Camera

- 5.0M CMOS High Resolution Camera
- 720P HD Video Up To 30fps
- True Color Restore
- Auto White Balance
- Auto/Manual Exposure
- Capturing Software, Free Driver Included

www.cnoec.com

