



CE / Rohs Biological Microscope Monocular Microscopes 400x Magnification

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1136
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **high powered microscope ,
binocular biological microscope**



for more products please visit us on cnoec.com

Product Description

CE / Rohs Biological Microscope Monocular Microscopes 400x Magnifica

A11.1136 Student Microscope,40x-400x

- 40x~400x Magnification
- LED Light Source
- Full Metal Body for Steday Observation
- Coaxial Coarse & Fine Focusing Knob

| A11.1136 Student Biological Microscope,40x-400x | |
|---|---|
| Head | Monocular Head |
| Eyepiece | WF10x |
| Objective | Achromatic 4x,10x,40x Spring |
| Nosepiece | Triple Revolving Nosepiece |
| Stage | Double Layers Mechanical Working Stage |
| Condenser | Abbe Condenser N.A.1.25 |
| Focusing | Coaxial Coarse & Fine Focusing Knob |
| Illumination | LED Light Source, Brightness Adjustable |



A11.1136 Biological Microscope



www.cnoec.com



A11.1136 Biological Microscope

www.cnoec.com



A11.1136 Biological Microscope

www.cnoec.com



A11.1136 Biological Microscope



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China