



## CCD Camera , Digital Microscope Camera Microscope Accessories A59.4903

Our Product Introduction

for more products please visit us on cnoec.com

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.4203
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Highlight: microscope accessory, compound microscope parts



### More Images



### Product Description

#### CCD Camera , Digital Microscope Camera Microscope Accessories A59.4903

**A59.4903 WIFI Digital Camera, 5.0M,iPad/Android/Win**,High quality upon very competitive cost.

- 5.0M CMOS High Resolution Camera
- WiFi IEEE 802.11b/g/n
- Frame Rate Up to 40 fps
- Support Microsoft Windows, MAC OS, IOS, Android or Above
- Enable Your Microscope Upgrade to WIFI Digital Microscope Easily!Auto/Manual Exposure

Our Product

--Capturing Software, Free Driver Included

A59.4903	A59.4903-M	A59.4903-M05	A59.4903-T	
Product Size	Dia. 60mm * Height 70mm			
Sensor	5.0M True Color 1/2.5"CMOS			
Resolution	2592 x 1944 ( 5.0M)			
Frequency Range	2.4000--2.4835GHz			
Working Modes	A.USB mode B.WiFi mode			
Focus Mode	Real-time Auto-focus			
Exposure	Wide Range Auto/manual Exposure			
White Balance	Automatic White Balance			
Function	Video Recording/ Audio Recording And Ttransmission, Photo Snapshot Function			
Turn On	WiFi Devices Instant Monitoring From Smartphone/ Tablet			
Support Browser	Safari / Chrome / IE			
Frame Rate	JPG Output: 640x480@40fps,800x600@40fps,1280x960@40fps, 1280x720@40fps, 2592x1944@15fps			
Transmit	Bandwidth 150M Bps			
WiFi	IEEE 802.11b/g/n			
Power supply	Separable and rechargeable lithium battery(3.7V 1900mAh); DC 5V USB power supply by PC External DC 5V 1.5A adapter			
Operating System	Microsoft® Windows® XP / Vista / 7 / 8 (32 & 64 bit) MAC OS X (Wi-Fi Working Mode) IOS5/IOS6(Wi-Fi Working Mode) Andriod 2.2 or above			
Hardware Requirement	CPU: Intel Core2 1.6GHz or Above Memory:2GB or Above USB Port:2.0 DVD-ROM(Required When Installing Software)			
Connecting	C-Mount	C-Mount, Including Adapter Ring 23.2mm, 30mm,30.5mm	C-Mount, Including Adapter Ring 0.5x/23.2mm, 30mm,30.5mm	C-Mount, Including Universal Telescope Adapter



**A59.4903 WIFI Digital Camera, 5.0M,iPad/Android/Win**





www.cnoec.com



**A59.4903 WIFI Digital Camera, 5.0M, iPad/Android/Win**  
With 5.0M High Resolution Digital Camera + WIFI Transfer Mode



www.cnoec.com



**A34.4903 WIFI Microscope**, Consist of  
**A59.4903 WIFI Digital Camera + A21.0902 Mono Video Microscope**



**A59.4903 WIFI Digital Camera, 5.0M, iPad/Android/Win**  
Compact Design, Mini Size, With Standard C-Mount Excellantly Work With Your  
Microscope or Telescope!

**A59.4903 WIFI Digital Camera, 5.0M, iPad/Android/Win**

With 5.0M High Resolution Digital Camera,  
WIFI Transfer Mode And U Cam Software,  
This WIFI Digital Camera Can Help You To  
Preview & Capture Image On Your iPhone,  
iPad or Android Mobile Phone!


You Can Easily View And Share The  
Excellent World Together With Your



Friends Now!

Mini USB Port

www.cnoec.com



www.cnoec.com



Back Auto exposure Zoom out Zoom in Auto focus Snapshot Video Folder  
Decrease exposure Increase exposure



www.cnoec.com



A59.4903 WIFI Digital Camera, 5.0M,iPad/Android/Win

A59.4903 WIFI Digital Camera, 5.0M,iPad/Android/Win



A59.4903 WIFI Digital Camera, 5.0M,iPad/Android/Win  
With 5.0M High Resolution Digital Camera + WIFI Transfer Mode



Mini USB Port


www.cnoec.com

A59.4903 WIFI Digital Camera, 5.0M, iPad/Android/Win



www.cnoec.com

A34.4903 WIFI Microscope, Consist of  
A59.4903 WIFI Digital Camera + A21.0902 Mono Video Microscope



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China