



## CCD 5.0M 720P Digital Microscope Camera 1/2.5" CMOS High Resolution A59.4904

Our Product Introduction

for more products please visit us on cnoec.com

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.4904
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Thery: Digital Camera
- Resolution: 5.0M
- Sensor: Cmos
- White Balance: Auto
- Interface: USB
- Lens Mount: C-Mount
- Highlight: microscope accessory, parts of the microscope



### More Images



## Product Description

CCD 5.0M 720P Digital Microscope Camera 1/2.5" CMOS True Color High Resolution

**A59.4904 Digital Camera, 720p, C Mount + Eyepiece Type**, High quality under very competitive cost.

- 5.0M CMOS High Resolution Camera
- 720P HD Video Up To 30fps
- True Color Restore
- Auto White Balance
- Auto/Manual Exposure
- Capturing Software, Free Driver Included


	A59.4904	A59.4904-M	A59.4904-M05
Digital Camera	5.0M CMOS High Resolution Digital Camera		
Senser	1/2.5" CMOS, True Color		
Resolution	2592*1944 Max.		
Pixel Size	2.2x2.2μm		
Frame Rate	@640X 480: 30fps; @2592X1944: 15fps 720P HD Video Up To 30fps		
White Balance	Auto		
Exposure	Manual/Auto		
Image Format	JPG/BMP/PNG/GIF		
Video Format	AVI		
Software	Capturing Software, Free Driver		
System	Microsoft Windows XP/2003/Vista/Win7/Win8		
Lens Mount	Standard C-mount		
Connecting	C-Mount	C-Mount, Including Adapter Ring 23.2mm, 30mm,30.5mm	C-Mount, Including Adapter Ring 0.5x/23.2mm, 30mm,30.5mm



A59.4904 Digital Camera, 5.0M, USB2.0




**Standard C-Mount**  
With Simple Installation Design, It Can Be Used in Trinocular Photo Tube With C-Mount Of All Kinds Biological Microscopes Or Stereo Microscopes



**5.0M High Resolution CMOS Sensor**  
Adopt The New CMOS Ture 5.0M Sensor, It Greatly Reduce Image Noise.



**A59.4904-M USB Digital Camera, 5.0M,720p**  
With Dia. 23.2mm Ring For Biological Microscope Eyepeice Tube  
With Dia. 30 and 30.5mm Ring For Stereo Microscope Eyepeice Tube



0.5X Lens Dia.23.2mm Is Equipped With Rubber Ring To Make Adapter Ring Very Stable  
Avoid Detachment Troubles



**A59.4904 Microscope Digital Camera**  
5.0M CMOS High Resolution Camera  
720P HD Video Up To 30fps  
True Color Restore  
Auto White Balance  
Auto/Manual Exposure  
Capturing Software, Free Driver Included

www.cnoec.com



A59.4904 Digital Camera, 5.0M, USB2.0




A59.4904, C-Mount

Camera 23.2mm 30mm 30.5mm

Camera 0.5x Lens 30mm 30.5mm

A59.4904-M

A59.4904-M05



0.5X Lens Dia.23.2mm Is Equipped With Rubber Ring To Make Adapter Ring Very Stable Avoid Detachment Troubles

## A59.4904 Microscope Digital Camera

5.0M CMOS High Resolution Camera  
720P HD Video Up To 30fps  
True Color Restore  
Auto White Balance  
Auto/Manual Exposure  
Capturing Software, Free Driver Included



[www.cnoec.com](http://www.cnoec.com)



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China