




Opto Edu 7" Lcd Handheld Digital Microscope Biological 2.0m

Our Product Introduction

for more products please visit us on cnoec.com


Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.5102
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Function: Wireless Digital Microscope Handheld
- LCD Screen: "7" IPS Screen, With Touch Button To Take Photo & Video, One Key Switch Between Main Body Camera & Portable USB Camera"
- Head: Mono Straight Head
- Eyepiece: WF10x, WF16x
- TF Card Output: Support Max 32G, To Take Photo & Video
- USB Output: To Win/Mac, With Professional StrangeView Measuring Software
- Power Supply: DC Adapter, Can Be Powered By USB Power Bank For Outdoor Use
- Package: Color Box + Foam, 1 Pc/1 Box, With Software Disc, DC Adapter, Instruction Manual
- Nosepiece: Triple
- Objective: Achromatic 4x10x40x
- Working Stage: Mechanical Stage



More Images



7" LCD Biological Digital Microscope, 2.0M

Dual Optical System: 7" LCD Digital Camera + USB Portable Digital Microscope

Total Magnification 40x-640x Achromatic Optical System

Double Layer Mechanical Stage 115x125mm With Disc Condenser

Coarse & Fine Focusing Knob With Insert LED Light Source Brightness Adjustable

With 10 Prepared Slides, With Experiment Kit, Pack In Hand Carry Case+Carton

A33.5102 7" LCD Biological Digital Microscope, 2.0M	
LCD Screen	7" IPS Screen, With Touch Button To Take Photo & Video, 2.0M, 30fps
Head	Mono Straight Head
Eyepiece	WF10x, WF16x
Nosepiece	Triple
Objective	Achromatic 4x10x40x
Working Stage	Mechanical Stage
Focusing	Coarse & Fine Focusing
Illumination	Insert LED
TF Card Output	Support Max 32G, To Take Photo & Video
USB Output	To Win/Mac, With Professional StrangeView Measuring Software
Power Supply	DC Adapter, Can Be Powered By USB Power Bank For Outdoor Use
Package	Color Box + Foam, 1 pc/1 Box, With Software Disc, DC Adapter, Instruction Manual



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel: +8610 88696020, E-mail: sale@optoedu.com

A33.5102

7" LCD Biological + USB Portable Dual Lens Digital Microscope, 2.0M




Ergonomics

The best viewing angle range, sitting for a long time will not be tired and will not hurt your eyes

7" LCD Measurable

The screen can adjust the visual angle of 180°, real-time real-time measurement of cell size





Objective & Stage

Triple, Achromatic 4x10x40x, Mechanical Stage



Focusing

Coarse & Fine Focusing, Brass dovetail guide rail, high precision and durable



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China