

for more products please visit us on cnoec.com



OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A12.2601

Laboratory Biological Microscope



Seidentopf Head Inclined 30°
Rotatable 360° , Interpupillary Distance 48-75mm

Seidentopf Head

Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, High Eyepoint Eyepiece,
Which Is Convenient For People With Different Sitting Height; The Integrated Design Of The Eyepiece Lens Barrel,
Avoids The Loss And Contamination Of The Eyepiece And Enhances The Integrated Performance Of The Microscope



Exclusive And Innovative Low-power Dimming Objective

Quadruple, Backward

4 Sets Objective Available: Achromatic, Plan, Infinity E-Plan, Infinity Plan.



Mechanical Stage Size 145 x 140mm
Travel Range 76 x 52mm

Double Slide Holder Working Stage

Mechanical Stage Size: 145 mm x 140mm, Travel: 76mm x 52mm, Scale: 0.1mm



Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down



3W-LED Illumination Systems

Condenser & Light Source

Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down, With Filter Holder
3W-LED Illumination Systems



Coaxial Coarse & Fine Focusing



Holding Hole On Back



Wire Winding Device



Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use



Different Color Available To Customize Your Model!

Items	A12.2601 Biological Microscope Specification	-A	-AT	-B	-BT	-C	-CT	-D	-DT	Cata. No.
Head	Seidentopf Binocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm	●	○	●	○	●	○	●	○	A53.2621-B
	Seidentopf Trinocular Head Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, Light Split Fixed P20:E80	○	●	○	●	○	●	○	●	A53.2621-T
Eyepiece	WF10x/18mm	●●	●●							A51.2602-1018
	WF10x/20mm (High Eyepoint)			●●	●●	●●	●●	●●	●●	A51.2602-1020
Nosepiece	Quadruple, Backward	●	●	●	●	●	●	●	●	A54.2610-N734
	Quintuple, Backward	○	○	○	○	○	○	○	○	A54.2610-N735
Achromatic Objective	4x/0.1, W.D.=18.6mm	●	●							A52.2601-4
	10x/0.25, W.D.=6.5mm	●	●							A52.2601-10
	20x/0.40, W.D.=1.75mm	●	●							A52.2601-20
	40x/0.6(S), W.D.=0.47mm	●	●							A52.2601-40
	60x/0.80(S), W.D.=0.10mm	○	○							A52.2601-60
100x/1.25(S,O), W.D.=0.07mm	○	○							A52.2601-100	
Plan Objective	4x/0.1, W.D.=14.5mm			●	●					A52.2603-4
	10x/0.25, W.D.=5.65mm			●	●					A52.2603-10
	20x/0.40, W.D.=5.65mm			○	○					A52.2603-20
	40x/0.6(S), W.D.=0.85mm			●	●					A52.2603-40
	100x/1.25(S,O), W.D.=0.07mm			●	●					A52.2603-100
Infinity E-Plan Objective	EPLAN 4x/0.10, W.D.=12.1mm					●	●			A52.2605-4
	EPLAN 10x/0.25, W.D.=4.14mm					●	●			A52.2605-10
	EPLAN 20x/0.4					○	○			A52.2605-20
	EPLAN 40x/0.65(S), W.D.=0.58mm					●	●			A52.2605-40
Infinity Plan Objectives	EPLAN 100x/1.25(S, Oil), W.D.=0.19mm					●	●			A52.2605-100
	PLAN 2.5x/0.07, W.D.=8.47mm							○	○	A52.2606-2.5
	PLAN 4X/0.10, W.D.=12.1mm							●	●	A52.2606-4
	PLAN 10X/0.25, W.D.=4.64mm							●	●	A52.2606-10
	PLAN 20X/0.40(S), W.D.=2.41mm							○	○	A52.2606-20
	PLAN 40X/0.66(S), W.D.=0.65mm							●	●	A52.2606-40
	PLAN 60X/0.80(S), W.D.=0.33mm							○	○	A52.2606-60
	PLAN 100X/1.25(S, Oil), W.D.=0.12mm							●	●	A52.2606-100
	PLAN 100X/1.15(S, Water), W.D.=0.19mm							○	○	A52.2606-100W
Working Stage	Mechanical Stage Size: 145 mm×140mm, Travel: 76mm×52mm, Scale: 0.1mm, Double Slide Holder	●	●	●	●	●	●	●	●	A54.2601-S04
	Rackless (Integrated) Stage, Size 140×135mm, Travel Range 75×35mm, Double Slide Holder	○	○	○	○	○	○	○	○	A54.2601-S12
Focusing	Coaxial Coarse & Fine Focusing, Coarse Focusing Travel Range: 26mm, Fine Focusing Travel Range: 2μm/Division	●	●	●	●	●	●	●	●	
Condenser	Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down, With Filter Holder	●	●	●	●	●	●	●	●	A56.2614-C06b
Illumination	3W-LED Illumination Systems Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	●	●	●	●	●	●	●	●	A56.2613-3W
	12V/20W Halogen Lamp, Power Supply Wide Voltage 100V-240V	○	○	○	○	○	○	○	○	A56.2650-20W
	6V/30W Halogen Lamp, Power Supply Wide Voltage 100V-240V	○	○	○	○	○	○	○	○	A56.2650-30W
Field Diaphragm	Upgrade To Kohler Illumnation	○	○	○	○	○	○	○	○	A56.2615-S
Filter	Dia.45mm, Blue	○	○	○	○	○	○	○	○	A56.2616-45B
	Dia.45mm, Green	○	○	○	○	○	○	○	○	A56.2616-45G
	Dia.45mm, Yellow	○	○	○	○	○	○	○	○	A56.2616-45A
Dark Field	Dark Field Condenser, Dry, N.A.0.83-0.91, For 4x10x40x Objective	○	○	○	○	○	○	○	○	A5D.2610-S
Polarizing	Analyzer + Polarizer Set	○	○	○	○	○	○	○	○	A5P.2601-S
Phase Contrast	Infinity Plan Phase Contrast Objective 10x + Phase Contrast Slide 10x + Phase Contrast Condenser + Telescope 11x + Green Filter	○	○	○	○	○	○	○	○	A5C.2601-10
	Infinity Plan Phase Contrast Objective 20x + Phase Contrast Slide 20x + Phase Contrast Condenser + Telescope 11x + Green Filter	○	○	○	○	○	○	○	○	A5C.2601-20
	Infinity Plan Phase Contrast Objective 40x + Phase Contrast Slide 40x + Phase Contrast Condenser + Telescope 11x + Green Filter	○	○	○	○	○	○	○	○	A5C.2601-40
	Infinity Plan Phase Contrast Objective 100x + Phase Contrast Slide 100x + Phase Contrast Condenser + Telescope 11x + Green Filter	○	○	○	○	○	○	○	○	A5C.2601-100
	Turret Phase Contrast Attachment Infinity Plan Phase Contrast Objective 10x20x40x100x, Turret Phase Contrast Condenser 5-Holes, Telescope 11x, Green Filter	○	○	○	○	○	○	○	○	A5C.2603-S
	Upgrade to A16.2601-B2/T2 Fluorescence Microscope With: B,G China Fluorescent Set 100W Mercury Lamp House 100W Mercury Power Box, Digital Disp. 6 Holes Disc Fluorescent Unit Ultraviolet Protection Barrier Centering Objective China 100W Mercury Lamp 2 pcs						○	○		A5F.2601-2




A52.2601
Achromatic Objective



A55.2610
Digital Eyepiece
Camera Adapter



A51.2603
Plan Objective



A56.2615-S
Field Diaphragm



A55.2601-10
C-Mount 1.0x,
Focus Adjustable



A5P.2601-S
Analyzer + Polarizer Set



A56.2601
Reflecting Mirror & Holder



A55.2601-75
C-Mount 0.75x,
Focus Adjustable



A5D.2610-S
Dark Field Condenser, Dry



A55.2601-05
C-Mount 0.5x,
Focus Adjustable



A54.2601-S12
Rackless (Integrated) Stage



A56.2616
Filter : Blue, Green, Yellow



A5C.2601
Phase Contrast Objective + Side + Condenser + Telescope + Green Filter



A5C.2603-S
Turret Phase Contrast Attachment

Product Multi-angle






Product Upgrade



A16.2601 Fluorescence Microscope



A16.2601-L LED Fluorescence Microscope



A13.2605 Metallurgical Microscope(Reflect)



A15.2603 Polarizing Microscope(Transmit)



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China