



Compensation Binocular Optical Microscope Halogen Illumination Microscopes

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1102
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **binocular compound microscope, wide field microscope**

for more products please visit us on cnoec.com

Product Description

Compensation Binocular Optical Microscope Halogen Illumination Microscopes

A12.1102 Biological Microscope		-A	-B	-C
Head	Monocular Head,45° inclined,360° Rotatable	●	-	-
	Compensation Free Binocular,30°inclined,55mm-75mm	-	●	-
	Compensation Free Trinocular,30°inclined,55mm-75mm	-	-	●
Eyepiece	WF10X/18mm	●	●	●
Objective	Achromatic 4X	●	●	●
	Achromatic 10X	●	●	●
	Achromatic 40X, Spring	●	●	●
	Achromatic 100X, Spring,Oil	●	●	●
Nosepiece	Backward Quadruple Revolving Nosepiece	●	●	●
Stage	Double Layer Mechanical Stage,Size140x132mm,Moving Range:70mmx30mm	●	●	●
Focusing	Coaxial Coarse and Fine Focusing Scale Value 0.002mm	●	●	●
Condenser	ABBE N.A.1.25Condenser with Iris Diaphragm & Filter	●	●	●
Light Source	Halogen Bulb6V/20W or 6V/15W,AC 220V/110V, Adjustable Brightness	●	●	●
Optional Accessories				
Eyepiece	WF16X			
	WF20X			
	P16X			
Objective	Achromatic 20X			
	Achromatic 60X, Spring			
	Semi-plan Objective 4X			
	Semi-plan Objective 10X			
	Semi-plan Objective 20X			
	Semi-plan Objective 40X(S)			
	Semi-plan Objective 60X(S)			
	Semi-plan Objective 100X(S)			
	Plan Objective 4X			
	Plan Objective 10X			
	Plan Objective 20X			
	Plan Objective 40X(S)			
	Plan Objective 60X(S)			
Plan Objective 100X(S)				
Dark Field	Darkfield			
Photo Adapter	Photo Adapter			
CCD Adapter	0.5x			
	0.57x			
	0.75x			



A12.1102-A

A12.1102-B



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China