



8' LCD Digital Eyepiece Camera Microscope Accessories A59.5601

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.5601
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: microscope accessory, compound microscope parts



More Images



Product Description

8" LCD Digital Eyepiece Camera Microscope Accessories A59.5601

A59.5601 8" LCD Pad Digital Camera,

- 8" LCD Touch Screen, 1024*768
- 2.0M High Resloution Digital Camera Built-in
- HDMI, USB2.0, Mini USB 2.0 Output
- Installed Measuring Software, Capture Picture & Video
- Input AC 110~220V Wide Voltage, Output 5V

A59.5601 8" LCD Pad Digital Camera	
LCD	8" High Resolution LCD, Color, Touch Screen, 1024*768
Camera	2.0M Digital Camera Built-in
FPS	15fps @1024*768
System	Android 4.1.1
CPU	Dual Core 1.5GHz
Memory	1GB
Storage	Internal 5GB, Extended By Mini SD Card
Output	HDMI, USB 2.0, Mini USB 2.0
Power	Input AC 110~220V Wide Voltage, Output 5V



A59.5601 8" LCD Pad Digital Camera, 2.0M



Measuring Software, Capture Picture / Video

www.cnoec.com



A59.5601 8" LCD Pad Digital Camera, 2.0M, Back View



U-TV C-Mount To Use On Microscope Body

www.cnoec.com



A59.5601 8" LCD Pad Digital Camera, 2.0M
With U-TV C-Mount To Use On Microscope Body



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China