



## HDMI 2.0M Digital Microscope Camera 1/2.8" CMOS Color Progressive Exposure

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.4209
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Highlight: microscope accessory, parts of the microscope



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

HDMI 2.0M Digital Microscope Camera 1/2.8" CMOS Color Progressive Exposure

<b>A59.4209 HDMI Digital Camera</b>	
Senser	1/2.8" CMOS , Color, 2.0M
Pixel Size	2.8×2.8μm
Resolution	1080P@60Hz
Output Digi	8 Bits
Exposure	Progressive Exposure
SNR	>52dB
FPS Max	30fps
Output	HDMI
Power Supply	DC 12V
Max Power	3W
Working Temperature	0°C to +50°C
Lens Mount	Standard C/CS Mount
Dimension	51mm×52mm×65mm
External Control	Remote Control
Basic Function	--Frezze --Auto White Balance --Auto Exporsure --Brightness Adjustable --8 Crosshairs On Position Display,Adjustable Colors -- Black and Wihte video change



A59.4209 HDMI Digital Camera



A59.4209 HDMI Digital Camera Back View



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

