



Heating Stage Portable Digital Microscope 40x H10x Eyepiece

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.5101
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Function: Handheld Portable Digital Microscope
- Screen: 7" LCD
- Magnification: 40x-640x
- Eyepiece: H10x, H16x
- Head: Dual Viewing Head, 45°Including 360°Rotable
- Nosepiece: Triple
- Objective: Achromatic 4x, 10x, 40x(s)
- Focusing: Coarse Distance: 20mm
- Stage: Heated Stage 37°C-38°C, Size 110mmx120mm, Moving 60mmx30mm
- Condenser: Rotating Diaphragm, $\Phi 3.2$ $\Phi 4$ $\Phi 5$ $\Phi 8$ $\Phi 15$
- Light Source: LED, 12V, Brightness Adjustable
- Highlight: heating stage portable digital microscope, 40x portable digital microscope, h10x eyepiece hand held digital microscope



More Images



Product Description

7" LCD Digital Heating Stage Biological Microscope

7" LCD Monitor + Digital Camera All-in-One For Easy Use

Heating Stage Temperature Adjustable Range 30°~50° For Blood, Sperm View

Total Optical Magnification 40x~640x, Digital Magnification 200x~2000x On Screen

Coaxial Coarse & Fine Focusing Knob, Focus Range 20mm

Transmit Light Source LED Brightness Powered By 12V DC Adapter

A33.5101 7" LCD Heating Stage Biological Microscope	
Screen	7" LCD
Magnification	40x-640x
Eyepiece	H10x, H16x
Head	Dual Viewing Head, 45°Including 360°Rotable
Nosepiece	Triple
Objective	Achroamtic 4x, 10x, 40x(s)
Focusing	Coarse Distance: 20mm
Stage	Heated Stage 37°C-38°C, Size 110mmx120mm, Moving 60mmx30mm
Condenser	Rotating Diaphragm, Φ3.2,Φ4,Φ5,Φ8,Φ15
Light Source	LED, 12V, Brightness Adjustable
Packing	Aluminium Box + Strong Master Carton

OPTO-EDU (BEIJING) CO., LTD.



A33.5101


7" LCD Heating Stage Biological Microscope

7" LCD, Dual
Viewing Head

Eyepiece
H10x, H16x

Triple Nosepiece,
Achromatic
4x10x40x

Heated Working Stage
37~38 °C




Coarse & Fine
Focusing Knob



Product Details



Achromatic 4x, 10x, 40x(s)



Heated Stage 37°C-38°C
Size 110x120mm, Moving Range 60x30mm



Temperature Control And Display Screen To
Ensure The Activity Of The Observed Sample
And Make The Observation Result More Accurate



Aluminum Box + Strong Master Carton



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China