



3.5" LCD Digital Microscope Accessories Eyepiece Camera , A59.1501

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.1501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **compound microscope parts,
parts of the microscope**



More Images



for more products please visit us on cnoec.com

Product Description

3.5" LCD Digital Microscope Accessories Eyepiece Camera , A59.1501

A59.1501 3.5" LCD Digital Eyepiece Camera

- 3.5" LCD Screen, Resolution Up To 12M(4023*3024)
- 2.0M CMOS Sensor, 1600*1200
- Output : SD Card Slot , USB 2.0, AV Output
- Time Stamp Function
- Standard Eyepiece Tube Size Dia. 23.2mm
- Multi Language Version, Default English, 9 Kinds Languages Version Built-in


A59.1501 3.5" LCD Digital Eyepiece Camera	
Screen	3.5" LCD
Resolution	2.0M CMOS Sensor, 1600*1200
LCD Resolution	12M(4023*3024) 8M(3264*2448) 5M(2560*1290) 3M(2048*1536) 2M(1600*1200) 1.3(1280*1024) VGA(640*480)
Memory	Internal Memory 128MB
Output	SD Card Slot
	USB 2.0 Output
	TV Output, NTSC/PAL
Quality	Fine, Standard, Economy
EV	Can Select EV from -2.0----+2.0
White Balance	Default Auto White Balance
ISO	Default Auto ISO
Color	Standard, Vivid, Sepia, Monochrome
Saturation	Normal, High, Low
Sharpness	Sharp, Normal Soft
Time Stamp	Off
	Date Only (Can Show The Date On The Picture Or Video)
	Date & Time (Can Show The Date & Time On The Picture Or Video)
Language	Default English, 9 Kinds Languages Version
Connect	Standard Eyepiece Tube Size Dia. 23.2mm



www.cnoec.com



A59.1501 3.5" LCD Digital Eyepiece Camera
Ideal Tools to Use With Biological Microscope In Eyepiece Tube,
Show Micro Images Clear For Share And Teaching



www.cnoec.com

A59.1501 3.5" LCD Digital Camera
SD Card Slot
USB 2.0 Output
AV Output



www.cnoec.com


A59.1501 3.5" LCD Digital Eyepiece Camera
Standard Eyepiece Tube Size Dia. 23.2mm For All
Biological Microscope



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China