



60 Microscope LED Ring Light Microscope Accessories Adjustable with UV Light

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.5202
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **microscope accessory, compound microscope parts**



for more products please visit us on cnoec.com

Product Description

60 Microscope LED Ring Light Microscope Accessories Adjustable with UV Light

A56.5202 Color LED Ring Light, 60 LEDs

- 60 LEDs Ring Light
- Four Colors Available : Blue, Red, Green, UV Light
- Brightness Adjustable 9~100%
- Outer Diameter 95mm, Inner Diameter 60mm
- Wide Range Voltage 110-265V

A56.5202 Color LED Ring Light	B	G	R	UV
Color	Blue	Green	Red	UV Light
Resolution	2.0M CMOS Sensor, 1600*1200			
LED	60pcs			
Angle	45°			
Outer Diameter	95mm			
Inner Diameter	60mm			
Height	30mm			
Power	≤5W			
Input Voltage	Wide Range Voltage 110-265V			
Output Voltage	DC13V			
Brightness Adjustment	0-100%			
Focal Length	100mm			
Working Distance	Effective Working Distance 50-200mm			



A56.5202 Color LED Ring Light,
Four Colors Available: Red, Green, Blue, UV Light



A56.5202-B

A56.5202-R

A56.2020-G



Front View

Back View



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China