



Trinocular Head Polarizing Light Microscope With Transmit Light CE A15.1303

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.1303
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **student compound microscope ,
highest magnification microscope**

for more products please visit us on cnoec.com

Product Description

Trinocular Head Polarizing Light Microscope With Transmit Light CE A15.1303

A15.1303 Transmit Polarizing Microscope, Newly designed with comfortable operating with ergonomic structure, with better illumination, convenient replacement of attachment with sliding-in center adjustable condenser.

A15.1303 Polarizing Microscope	
Head	Sliding Trinocular Head, 45° Inclined, 360° Rotatable, Interpupillary Distance 55~75mm
Eyepiece	WF10x/18mm
	Wide Field Plane-Scope Eyepiece WF10X/18 with scale of Crosshair and 0.1mm
Nosepiece	Backward, Quaduple Nosepiece
Objective	45mm Strain Free Achromatic 4x, 10x, 40x(S), 63x
Working Stage	Revolving Round Stage: Ø160mm, 360° Rotatable and Graduated 10 Increments, Minimum Resolution 6' When Using Vernier Scale Center Adjustable;
Condenser	Abbe N.A. 1.25, With Dia. 2-Dia 30mm Iris Diaphragm and Ø32 Filter
Focusing	Coaxial Coarse And Fine Focusing Adjustable, Mechanism: 35mm, Precision: 0.002mm;
Light Source	Adjustable Brightness Halogen Lamp 6V/20W (or 3.5V/1W LED Lamp) and Kohler Illumination System;
Polarizing	Rotatable Analyzer with Graduation 0°-90°, Sliding out of Optical Path;
	Bertrand Lens, Sliding out of Optical Path;
	Optical Compensator: λ Slip (first class red), 1/4 λ Slip, Quartz Wedge
Polarize:	Under the Stage, 0° Adjustable 360° Rotatable
Package	1 set/ Carton 42X28X57cm
Gross Weight	13kg






Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China