



OPTO-EDU A13.2602-B Inverted Metallurgical Microscope, Binocular, Infinity Plan

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.2602-B
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **digital trinocular microscope, usb digital microscope**



for more products please visit us on cnoec.com

Inverted Metallurgical Microscope, Binocular, Infinity Plan

Classical Best Compact & Economical Inverted Metallurgical Microscope!

Quadruple Nose With Metallurgical LWD Infinity Plan Objective

Mechanical Working Stage 180x180mm Moving Range 30x30mm

3W LED Standard Light Source, 12V20W Halogen Optional

Items	A13.2602 Inverted Metallurgical Microscope Specification	-A	-B	Cata. No.
Head	Monocular Head Inclined 45°, Rotatable 360°	●		A53.2631-M
	Seidentopf Binocular Head Inclined 45°, Rotatable 360°, Interpupillary Distance 48-75mm		●	A53.2631-B
Eyepiece	WF10x/18mm	●	●●	A51.2602-1018
	WF10x/18mm With Micrometer X-0.1mm*180 Div / Y-Single Line, Diopter Adjustable	○	○	A51.2605-1018D
Nosepiece	Quadruple	●	●	A54.2610-N704R
Metallurgical LWD Infinity Plan Objectives	LPL 5x/0.13, W.D.=16.04mm	●	●	A5M.2633-5
	LPL10x/0.25, W.D.=18.48mm	●	●	A5M.263310
	LPL 20x/0.40, W.D.=8.35mm	○	○	A5M.2633-20
	LPL 50x/0.70, W.D.=1.95mm	●	●	A5M.2633-50
	LPL 80x/0.80, W.D.=0.85mm	○	○	A5M.2633-80
	LPL 100x/0.90(Dry), W.D.=1.00mm	○	○	A5M.2633-100
Working Stage	Stage Size 180x180mm, Travel Range 30x30mm, Round Plate Size Dia.110mm	●	●	A54.2601-MP100
	Round Plate With Oval Hole	○	○	A54.2601-MD400
Light Source	3W LED Lamp, Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	●	●	A56.2612-3W
	12V 20W Halogen Lamp, Wide Voltage: 100V-240V	○	○	A56.2650-20W
Filter	Blue, Dia.32mm	●	●	A56.2616-32B
	Green, Dia.32mm	●	●	A56.2616-32G
	Grey, Dia.32mm	●	●	A56.2616-32N
Micrometer	Division 0.01mm	○	○	A58.2001-901
Package	Carton Size 437*268*522mm, 1 Unit/Ctn			
Weight	Gross Weight 8.0kgs, Net Weight 7.0kgs			
Note: "●" In Table Is Standard outfits, "○" Is Optional Accessories.				

Eyepiece Micrometer + Stage Micrometer For Measuring Function Optional

Items	A13.2602 Inverted Metallurgical Microscope Specification	-A	-B	Cata. No.
Head	Monocular Head Inclined 45°, Rotatable 360°	●		A53.2631-M
	Seidentopf Binocular Head Inclined 45°, Rotatable 360°, Interpupillary Distance 48-75mm		●	A53.2631-B
Eyepiece	WF10x/18mm	●	●●	A51.2602-1018
	WF10x/18mm With Micrometer X-0.1mm*180 Div / Y-Single Line, Diopter Adjustable	○	○	A51.2605-1018D
Nosepiece	Quadruple	●	●	A54.2610-N704R
Metallurgical LWD Infinity Plan Objectives	LPL 5x/0.13, W.D.=16.04mm	●	●	A5M.2633-5
	LPL10x/0.25, W.D.=18.48mm	●	●	A5M.263310
	LPL 20x/0.40, W.D.=8.35mm	○	○	A5M.2633-20
	LPL 50x/0.70, W.D.=1.95mm	●	●	A5M.2633-50
	LPL 80x/0.80, W.D.=0.85mm	○	○	A5M.2633-80
	LPL 100x/0.90(Dry), W.D.=1.00mm	○	○	A5M.2633-100
Working Stage	Stage Size 180x180mm, Travel Range 30x30mm, Round Plate Size Dia.110mm	●	●	A54.2601-MP100
	Round Plate With Oval Hole	○	○	A54.2601-MD400

Light Source	3W LED Lamp, Powered By Safety Low Voltage Charger, Input Wide Voltage 100V-240V, Output 5V1A, Support Power Bank Supply For Outdoor Use	•	•	A56.2612-3W
	12V 20W Halogen Lamp, Wide Voltage: 100V-240V	○	○	A56.2650-20W
Filter	Blue, Dia.32mm	•	•	A56.2616-32B
	Green, Dia.32mm	•	•	A56.2616-32G
	Grey, Dia.32mm	•	•	A56.2616-32N
Micrometer	Division 0.01mm	○	○	A58.2001-901



OPTO-EDU (BEIJING) CO., LTD.

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085

A13.2602 Inverted Metallurgical Microscope




Product Models & Size(mm)




A13.2602-A



A13.2602-B



A13.2602-A General Layout Diagram



A13.2602 Inverted Metallurgical Microscope Size(mm)

Product Accessories



A51.2602-1018
WF10x/18mm




A56.2616
Filter: Blue, Green, Grey



A54.2610-N704R
Quadruple




A5M.2633
Metallurgical LWD Infinity Plan Objectives



A54.2601-MD400
Round Plate With Oval Hole



A54.2601-MP100
Working Stage Size 180x180mm



A58.2001-901
Division 0.01mm Micrometer



A58.2002
Eyepiece Micrometer For Measuring Function



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China