



Extra Wide Field Stereo Optical Microscope With Big Base A22.0302 Ce Listed

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.0302
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Product Name: Stereo Microscope
- Application: Electrical Industry And So On
- Objective: 2x, 4x, 6x, 8x, 10x
- Base: 400×300×20mm
- Vertical Arm Height: 500mm
- Illumination: Halogen Lamp & LED Light
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



More Images



Product Description

Extra Wide Field Stereo Optical Microscope with Big Base A22.0302 With CE

- Original creation & designing in combination with visual project
- Easy changable objective 2x, 4x, 6x, 8x, 10x, humanisic operatiing system
- Suitable for long time working, comfortable operating and viewing improve efficiency
- Super long working distance, super large field of viewing & depth of field.
- Used in electrical industry, precision project, medical diagnosis and so on.

| A22.0302 Stereo Microscope | | | |
|----------------------------|----|----|----|
| Magnification | 2x | 4x | 6x |

Our Product

| | | | | |
|---|---|-----------------------|--------------------|-------------|
| Objective | Working Distance (mm) | 208 | 98 | 80 |
| | Viewing Field (mm) | 68 | 34 | 22.7 |
| Base Plate Size | 400x300x20mm | | | |
| Vertical arm height | 500mm | | | |
| Travel | 200mm | | | |
| Illumination | AC 110V/220V, 12V 20W halogen lamp with cup (2 pcs), Outside LED light adjuster | | | |
| A12.0205 Laboratory Biological Microscope Optional Accessories | | | | |
| Objective | Magnification | Working Distance (mm) | Viewing Field (mm) | Item No. |
| | 8x | 58 | 17 | A52.0380-08 |
| | 10x | 46 | 13.6 | A52.0380-10 |



E22.0302



E22.0302

www.cnoec.com



E22.0302



Large View Field Glass



Objective



Observing a Screw



Image Effect

www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China