



Portable / Mini Metallurgical Optical Microscope 100 - 400x A13.2501

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.2501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: reflected light microscope, bright field microscope



More Images



Product Description

Portable / Mini Metallurgical Optical Microscope 100 - 500x A13.2501
Small & Simple, with bright white LED light built-in, and high quality image!

A13.2501 Portable Metallurgical Microscope		
Total Magnification	100-400x	
Head	Monocular	
Eyepiece	WF10X	
LWD Plan Objective (no cover glass)	PLL 10X/0.25 WD:7.3mm PLL 40X/0.65 WD:0.5mm	
Mechanical Tube Length	160mm	
Focusing	Coaxial coarse/fine focus system, range:36mm	
Light Source	Coaxial illumination, build-in LED rechargeable, brightness adjustable	
Packing	Portable aluminium case	
A13.2501 Portable Metallurgical Microscope Optional Accessories		
Adapter	For Canon SLR Digital Camera	A55.2502
Stand	Magnet Stand	A54.2501
LWD Plan Objective (no cover glass)	PLL 20x	A5M.0234-20
	PLL 60x	A5M.0234-60
	PLL 80x	A5M.0234-80

Our Product Introduction

for more products please visit us on [cnoec.com](http://www.cnoec.com)



A13.2501





LWD Plan Objective



Focusing



Portable Aluminum Case



Optional Part: A54.2501 Magnet Stand



A55.2502 Canon SLR Camera Adapter



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627

sale@optoedu.com

cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China