



40X - 1000X Infinity Inverted Halogen Phase Contrast Microscope A19.2703

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A19.2703
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **phase contrast light microscope,
phase contrast microscopes**

Product Description

Infinity Achromatic Halogen Phase Contrast Microscope 40X - 1000X A19.2703

Specification	A19.2703 Inverted Phase Contrast Microscope
Optical System	UCIS Infinity independent achromatic optical system
Magnification	40x-1000x
Head	Seidentopf trinocular head 45° inclined, interpupillary distance 52-75mm
Eyepiece	WF10X/20mm, high eye-point up to 21mm, ±5 Diopter adjustment
Objective	Infinity Achromatic objective 4x
	Infinity Achromatic phase-contrast objective 10x
	Infinity Achromatic phase-contrast objective 40x(S)
	Infinity Achromatic phase-contrast objective 100(S/Oil)
Focus System	Coaxial coarse and fine focus mechanism with markings. Fine focus sensitivity 0.001mm
Nosepiece	Quintuple nosepiece inward facing with positive click stops
Condenser	N.A. 0.30 Abbe condenser with iris diaphragm, with phase-contrast slide, long work distance 72mm
Stage	Right hand controls, brightness adjustable, low position, mechanical stage, 240mm x 260mm, movement range 135mm x 85mm
Adapter	0.5x C-mount
Light source	6V 30 W halogen illumination



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

