



Halogen Bulb Metallurgical Optical Microscope Infinity Objective A13.1109

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.1109
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **reflected light microscope, halogen lamp microscope**



More Images



for more products please visit us on cnoec.com

Product Description

Metallurgical Optical Microscope Infinity Objective Halogen Bulb A13.1109

| A13.1109 Metallurgical Microscope Specification | | -A | -B | -C | -D |
|---|---|----|----|----|----|
| Head | Compensation Free Binocular,30°inclined,55mm-75mm | ● | - | ● | - |
| | Compensation Free Trinocular,30°inclined,55mm-75mm | - | ● | - | ● |
| Eyepiece | Wide field WF10X/20mm | ● | ● | - | - |
| | Wide field WF10X/22mm | - | - | ● | ● |
| | Wide field WF10X/20mm with Crosshair 0.1mm | - | ● | - | ● |
| 195 Metallurgical Objective | 4X/0.1 W.D. 25mm | ● | ● | - | - |
| | 10X/0.25 W.D. 11mm | ● | ● | - | - |
| | 20X/0.4 W.D. 9mm | ● | ● | - | - |
| Infinity Metallurgical Objective | 40X/0.6 W.D. 3.8mm | ● | ● | - | - |
| | 4X/0.1 W.D. 25mm | - | - | ● | ● |
| | 10X/0.25 W.D. 11mm | - | - | ● | ● |
| | 20X/0.4 W.D. 8mm | - | - | ● | ● |
| Flashboard Filter | Blue Filter | ● | ● | ● | ● |
| | Green Filter | ● | ● | ● | ● |
| | Yellow Filter | ● | ● | ● | ● |
| Focus System | Coaxial Coarse and Fine Focusing with Rack and Pinion Mechanism; Fine Focusing Scale Value 0.002mm. | ● | ● | ● | ● |
| Condenser | ABBE N.A.1.25 Condenser with Iris Diaphragm & Filter | ● | ● | ● | ● |
| Nosepiece | Quadruple | ● | ● | ● | ● |
| Stage | Double Layer Mechanical, Size:180mmX150mm,Moving Range:75mmX50mm) | ● | ● | ● | ● |
| Illumination | EPI-Kohler Illumination with Apertrue Iris Diaphragm and Field Iris Diaphragm | ● | ● | ● | ● |
| | Transmission Light,Adjustble Brightness | ● | ● | ● | ● |
| Light Source | 12V/30W,AC 85V-230V Halogen Bulb | ● | ● | ● | ● |
| Polarizing Device | Analyzer rotatable 360°, Polarize &Analyzer can be Mmoved in/out of the Optical Path | ● | ● | ● | ● |
| Checking Tool | 0.01mm Micrometer | ● | ● | ● | ● |
| Optional Accessories | | | | | |
| Objective | 50X/0.75 W.D. 1.9mm | | | | |
| | 80X/0.9 W.D. 0.9mm | | | | |
| Photo Adapter | Photo Adapter | | | | |



A13.1109



Quadruple Nosepiece



Compensation Free Binocular
55mm-75mm



ABBE N.A. 1.25 Condenser with Iris
Diaphragm & Filter



Fine Focusing Scale Value 0.002mm



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China

