



Fluorescent Microscope Trinocular Head 40 - 1000X with CE A16.1103

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A16.1103
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **epi-fluorescence microscope, inverted fluorescence microscope**



More Images



for more products please visit us on cnoec.com

Product Description

Fluorescent Microscope Trinocular Head 40 - 1000X with CE A16.1103

A16.1103 Fluorescence Microscope Specification		A	B	C	D	
Head	Compensation Free Binocular,30°inclined	-	-	•	-	
	Compensation Free Trinocular,30°inclined	-	-	-	•	
	Sliding Binocular,45°inclined	•	-	-	-	
	Sliding Trinocular,45°inclined	-	•	-	-	
Eyepiece	Wide field WF10X/18mm	•	•	-	-	
	Wide field WF10X/22mm	-	-	•	•	
Objective	195 Fluorescent Objective	4X	•	•	-	
		10X	•	•	-	
		40X(S)	•	•	-	
		100(S)Oil	•	•	-	
	Infinity Fluorescent Objective	4X	-	-	•	•
		10X	-	-	•	•
		40X(S)	-	-	•	•
		100(S)oil	-	-	•	•
Focus System	Coaxial Coarse and Fine Focusing with Rack and Pinion Mechanism; Fine Focusing Scale Value 0.002mm.		•	•	•	•
Condenser	N.A. 1.25 Abbe Condenser with Iris Diaphragm&Filter		•	•	•	•
Stage	Double Layers Mechanical Stage,Stage Size:180x150,Moving Range:75x50		•	•	•	•
Illumination	Epi- Illumination		•	•	•	•
	Transmission Illumination		•	•	•	•
Light Source	Upper:Mercury Lamp 220V(110V)/100W Pointer or Digital Display Power Supply Box		•	•	•	•
	Halogen Bulb:12V/30W		•	•	•	•
			•	•	•	•
Fluoroscope Box	B Wave Band		•	•	•	•
	G Wave Band		•	•	•	•
Optional Accessories						
Objective	20X					
	60X(S)					
Photography Adapter	Photography Adapter					
Fluoroscope Box	U Wave Band					
	V Wave Band					



Trinocular head



Nosepiece and Stage



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China