



A23.1001-T LED Zoom Stereo Microscope With Digital Slr Camera

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1001-T
- Minimum Order Quantity: 1pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: Western Union, T/T Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Optical System: Parallel Optical System, Module Design To Add Function Easily
- Head: Tilting Head 0°-35° Viewing Angle Adjustable For User
- Eyepiece: Extra Wide Field Eyepiece EWF10x/22mm
- Auxiliary Lens: Plan Achromatic Lens 1.0x
- Working Distance: 78mm
- Focusing: Coaxial Coarse & Fine Focusing Knob, Fine Division 0.015mm, Coarse Stroke 38.4mm/Rotation, Fine Stroke 1.5mm/Rotation, Moving Range 105mm, Focusing Mount Inner Dia. 76mm
- Illumination: Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours
- Zoom Body: 0.8~5.0x, Zoom Ratio 1:6 0.8~6.4x, Zoom Ratio 1:8 0.8~8.0x, Zoom Ratio 1:10
- Highlight: **led type zoom stereo microscope,**

for more products please visit us on cnoec.com

Product Description

Product Specification

A23.1001-T Zoom Stereo Microscope, With Tilting head		A23.1001-T06	A23.1001-T08	A23.1001-T10
Optical System	Parallel Optical System, Module Design to Add Function Easily	●	●	●
Head	Tilting Head 0°-35° Viewing Angle Adjustable for User	●	●	●
Eyepiece	Extra Wide Field Eyepiece EWF10x/22mm	●	●	●
Zoom Body	0.8~5.0x, Zoom Ratio 1:6	●	-	-
	0.8~6.4x, Zoom Ratio 1:8	-	●	-
	0.8~8.0x, Zoom Ratio 1:10	-	-	●
Auxiliary Lens	Plan Achromatic Lens 1.0x	●	●	●
Working Distance	78mm	●	●	●
Focusing	Coarse Focusing Knob, Range 105mm	●	●	
	Coaxial Coarse & Fine Focusing Knob, Fine Division 0.015mm, Coarse Stroke 38.4mm/Rotation, Fine Stroke 1.5mm/Rotation, Moving Range 105mm, Focusing Mount Inner Dia. 76mm			●
Illumination	Transmit & Reflect LED Illumination, Brightness Adjustable, Working Life 6000 Hours			
A23.1001-T Zoom Stereo Microscope, With Tilting head		Cata. No.		
Eyepiece	WF15x/16mm	A51.1027-1516		
	WF20x/12mm	A51.1027-2012		
	WF30x	A51.1027-30		
Auxiliary Lens	Plan Apochromatic Lens 0.5x	A52.1023-05A		
	Plan Achromatic Lens 1.0x	A52.1023-10		
	Plan Achromatic Lens 2.0x	A52.1023-20		
	Achromatic Lens 0.3x	A52.1023-03		
	Achromatic Lens 0.5x	A52.1023-05		
Illumination	Coaxial Illumination Set:	A56.1021		
	--Cold Light Power Supplier,			
	--Dual Light Guide,			
	--Coaxial Episcopic Illuminator,			
	--Iris Diaphragm			
	Iris Diaphragm	A56.1021-2		
Working Stage	Temperature Control Stage	A54.1004		
	Movable Working Stage	A54.1007		
Stand	Universal Table Stand	A54.1011		
		A54.1012		
		A54.1013		
		A54.1014		
		A54.1019		
	Track Stand, Transmit & Reflect Light Source With Coarse & Fine Focusing Optional, Model -F, Like A54.1019-F			
		A54.1020		
		A54.1021		
Fluorescent	EPI-Fluorescent Attachment For Stereo Microscope	A5F.1020		
Dark Field	Dark Field Attachment For Stereo Microscope	A5D.1001		
Polarizing	Simple Polarizing Attachment For Stereo Microscope	A5P.1003		
Adapter	Photo Tube, Beamsplitter 1 Port	A55.1003-S1		
	Photo Tube, Beamsplitter 2 Port	A55.1003-S2		
	Photo Adapter, 2.0x, For Nikon or Canon Digital SLR Camera	A55.1020-2X		
	C-mount 0.5x	A55.1021-05		



A23.1001-F Flourescent Zoom Stereo Microscope



Product Details



Head And Eyepiece


A23.1001 With Binocular Head Inclined 20°
A23.1001-T With Tilting Binocular Head 0~35°



Product Details



Parallel Optical System Model Design to Add Function Easily



Product Details



A5F.1020 EPI-Fluorescent Attachment For Stereo Microscope



A23.1001
Zoom Stereo Microscope, Max 0.8~8x



A23.1001-F
Flourescent Zoom Stereo Microscope

Product Details




A23.1001-T Standard Set



A23.1001-T Standard Set with
A55.1003-S1 Photo Tube, Beam Splitter 1 Port

Product Dimension



A23.1001

A23.1001-T

A23.1001-F

Accessories

A56.1021-2 Coaxial Episcopic Iris Diaphragm
Coaxial Episcopic Illumination
Dual Light Guide Cold Light Power
A56.1021 Coaxial Illumination Set

A54.1004 Temperature Control Stage

A54.1007 Movable Working Stage

A5D.1001 Dark Field Attachment

A54.1014 Tilting Head 0° -35°

A5F.1020 EPI-Fluorescent Attachment

A59.1006 8.4" TFT Touching LCD&Connector

A5P.1003 Simple Polarizing Attachment

A23.1001-T

A23.1001-F

A36.1001

A54.1011 A54.1012 A54.1013 A54.1014

A54.1019 A54.1020 A54.1021

