



## A22.3670 OPTO EDU Stereo Optical Microscope 1x2x Or 1x3x Or 2X/4X Binocular

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.3670
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



### Product Specification

- Eyepiece: WF10x/23mm
- Head: 45°Binocular Head
- Objective: Turret Objective 1x2x Or 1x3x Or 2X/4X
- Magnification Range: 10X-120X
- Working Distance: 110mm
- Stand: Pole Stand, No Light Source
- Highlight: pole stand stereo optical microscope, binocular stereo optical microscope, pole stand opto edu microscope



### More Images



## Product Description

The stereoscopic microscope has a positive stereoscopic sense, which is easy to operate and dissect. The field of view is large in diameter, the working distance is long, and the depth of focus is large, which is convenient to observe the whole layer of the object being detected.

Binocular Head, 45° Inclined, 360° Rotatable

Wide Field, High-Eyepoint Eyepiece-WF10X23mm

Standard Turret Objective 2x/4x, Optional 1x/3x, 1x/2x

Total Magnification 20x/40x, Up To 10x-120x With Optional Eyepieces & Objectives

35 Kinds Stand Optional, Details See A54.3630

A22.3670-B1 Stereo Microscope Specification	
Head	45° Binocular head
Eyepiece	Wide Field, High-eyepoint Eyepiece-WF10X23mm
Objective	Turret Objective 1x2x or 1x3x or 2X/4X
Magnification Range	10X-120X
Working Distance	110mm
Stand	Pole Stand, No Light Source

## A22.3670 Details



Eyepiece WF10X/23mm Super Wide Field, High Eyepoint, Eyepiece Tube Diopter Adjustable



20x/40x, Magnification Up To 10x, 120x With Optional Eyepieces And 1X2X, 1X3X, 2X4X, Objectives

## A22.3670 Series Models



Binocular Head



35 Stands



A22.3670-B1



A22.3670-B2



A22.3670-B2L



A22.3670-B3



A22.3670-B3LS



A22.3670-B2LD



A22.3670-B4



A22.3670-B5



A22.3670-B6



A22.3670-B7



A22.3670-B8L



A22.3670-B8LD



A22.3670-B9L



A22.3670-B9LD



A22.3670-B8LS



A22.3670-L



A22.3670-2L



A22.3670-B9LS





A22.3670-STL6A



A22.3670-STL6B



A22.3670-STL6C

Objective Optics Parameter

Magnification	Eyepiece				Working Distance
	10X23	10X20	15X16	20X12	
1X	10X23mm	10X20mm	15X16mm	20X12mm	110mm
2X	20X11.5mm	20X10mm	30X8mm	40X6mm	
3X	30X7.7mm	30X6.7mm	45X5.3mm	60X4mm	
4X	40X5.8mm	40X5mm	60X4mm	80X3mm	



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China