



CMOS Color Microscope Accessories , 2.0M CCD Camera For Microscope A59.4209

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.4209
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



A59.4209 HDMI Digital Camera

Product Specification

- Highlight: **compound microscope parts, parts of the microscope**



A59.4209 HDMI Digital Camera

for more products please visit us on cnoec.com

More Images



A59.4209 HDMI Digital Camera



HDMI Digital

A59.4209 HDMI Digital Camera

Product Description

CMOS Color Microscope Accessories , 2.0M CCD Camera For Microscope A59.4209

Our Product Introduction

A59.4209 HDMI Digital Camera	
Senser	1/2.8" CMOS , Color, 2.0M
Pixel Size	2.8x2.8µm
Resolution	1080P@60Hz
Output Digi	8 Bits
Exposure	Progressive Exposure
SNR	>52dB
FPS Max	30fps
Output	HDMI
Power Supply	DC 12V
Max Power	3W
Working Temperature	0°C to +50°C
Lens Mount	Standard C/CS Mount
Storage	SD Card
Dimension	51mmx52mmx65mm
External Control	Remote Control
Basic Function	--Freeze --Auto White Balance --Auto Exposure --Brightness Adjustable ---8 Crosshairs On Position Display,Adjustable Colors -- Black and White video change

for more products please visit us on cnoec.com



A59.4209 HDMI Digital Camera






Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China