



Biological Microscope Binocular Microscopes halogen lamp A11.0206

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.0206
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **high powered microscope ,
biological compound microscope**

Product Description


Biological Microscope Binocular Microscopes halogen lamp A11.0206

This microscope is the conveniently operated optical instruments with excellent functions and new pattern it is very suitable for usage in medical test, teaching and microobservation in hospitals and schools.


--With wide field eyepieces and achromatic objectives, the field of view is widely and clear

--6V 20W halogen lamp with brightness control

Specification	A11.0206-A	A11.0206-B
Eyepiece	WF 10X(Φ18mm)	
Objective	Achromatic 4X/0.10	
	Achromatic 10X/0.25	
	Achromatic 40X/0.65(Spring)	
	Achromatic 100X/1.25(Spring, oil)	
Head	Sliding binocular head inclination 30°	
Focus System	Coaxial coarse/fine system focus, Minimum division of fine focusing: 4μm	Coaxial coarse/fine system focus, Minimum division of fine focusing: 4μm With tensional adjustable and up stop
Nosepiece	Quadruple backward ball bearing inner locating	
Stage	Double layer mechanical, size 125x115mm, moving range 70x30mm	
Condenser	Abbe N.A.1.25 Rack & pinion adjustment, with aperture diaphragm, filter seat	
Filter	Blue filter	
	Frosted filter	
Collector	For halogen lamp	
Light source	6V 20W halogen lamp with brightness control	
Optional Accessories		Item No.
Eyepiece	Wide field WF16x/11mm	A51.0203-16A
	Dividing 10x/18mm	A51.0205-10
Head	Monocular, inclination of 45°	A53.0201-4
	Sliding trinocular, inclination of 45°	A53.0204-4
	Seidentopf binocular head, inclination of 30°	A53.0203-2
Nosepiece	Triple nosepiece (Backward ball bearing inner locating)	A54.0201
Stage	Rotatable round stage Dia.120mm	A54.0206-1
	Fixed, size 120x120mm, with single lens and iris diaphragm	A54.0205-07
Filter	Green filter	A56.0209-2
	Yellow filter	A56.0209-3
Light Source	Fluorescent lamp 5W, 110V/220V	A56.0202-4
	Incandensent lamp	A56.0203-3
	LED Illumination, high brightness white	A56.0204-2
	Transformer for LED, Input 230V, 50/60Hz, Output 9V500mA	A56.0205-1
	Transformer for LED, Input 110V, 50/60Hz, Output 9V500mA	A56.0205-2
Polarizing Set	Reflecting mirror, plane-concave reflecting mirror with stand	A56.0206-5
	Polarizer 360° rotatable stage	A5P.0210-1
	Analyzer	A5P.0211-1



A11.0207-B



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China