



LED Microscope Accessories , 5X Mini Magnifier Handheld Microscope G11.5001

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: G11.5001
- Minimum Order Quantity: 1
- Packaging Details: Carton
- Supply Ability: 1000



Product Specification

- Highlight: **microscope accessory,
compound microscope parts**

Product Description

LED Microscope Accessories , 5X Mini Magnifier Handheld Microscope G11.5001

- 1.Magnification of lens: 5X
- 2.3 pcs LED
- 3.OEM

G11.5001 Magnifier, With LED, 5X	
Magnification	5x
Material	Aluminium Alloy /Plastic
Lens Dia.	7cm
Light source	3 pcs LED
Battery	5 hours battery working time
Weight	210g (include battery)
Size of Product	W92*D92*H52
Accessories	Special Bag, 1pc Lens Cloth, 1pc Battery, 1pc Manual



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627
 sale@optoedu.com
 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China