



Halogen Lamp Biological Microscope Binocular / Trinocular Microscopes A11.0302

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.0302
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **high powered microscope ,
binocular biological microscope**



More Images



for more products please visit us on cnoec.com

Product Description

Halogen Lamp Biological Microscope Binocular / Trinocular Microscopes A11.0302

- A11.0302 Student Microscope,
 --Achromatic Objectives
 --WF10x/18mm Wide Field Eyepiece
 --Monocular, Binocular, Trinocular Head Optional
 --Mechanical stage

Specification	A11.0302 Biological Microscope	M	D	B	T
Head	Monocular head, inclined 45°, rotatable 360°	•			
	Dual Viewing head, inclined 45°, rotatable 360°	•			
	Binocular head, inclined 30°, interpupillary distance 55-75mm			•	
	Trinocular head, inclined 30°, interpupillary distance 55-75mm				•
Eyepiece	WF10x/18mm				
Nosepiece	Quadruple nosepiece				
Objective	Achromatic Objective 4x/NA 0.1				
	Achromatic Objective 10x/NA 0.25				
	Achromatic Objective 40x(s)/NA 0.65				
	Achromatic Objective 100x(s,oil)/NA.1.25				
Stage	Mechanical stage size 140x120mm, moving range 70x40mm, minimum vernier division 0.1mm				
Focusing	Coaxial coarse & fine focusing knob, coarse focusing range 10mm, fine focusing division 0.002mm				
Condenser	N.A.1.25 Abbe condenser, up/down range 15mm				
Light Source	6V 20W Halogen lamp, brightness adjustable				
Optional Accessory		Item No.			
Eyepiece	WF16x	A51.0131-16			
	WF20x	A51.0131-20			
Objective	Achromatic 20x (S)	A52.0301-20			
	Achromatic 60x (S)	A52.0301-60			
	Plan 4x	A52.0303-04			
	Plan 10x	A52.0303-10			
	Plan 20x(S)	A52.0303-20			
	Plan 40x(S)	A52.0303-40			
	Plan 60x(S)	A52.0303-60			
	Plan 100x(O,S)	A52.0303-100			
Condenser	Dark field condenser, dry	A57.0310-01			
	Dark field condenser, immersion	A57.0310-02			



A11.0302-B



www.cnoec.com

A11.0302-T



www.cnoec.com

A11.0302-M

A11.0302-D



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

