



## Sliding / Seidentopf Binocular Biological Microscope with 6V 20W Halogen Lamp A11.0303

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.0303
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Highlight: **binocular biological microscope, biological compound microscope**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Sliding / Seidentopf Binocular Biological Microscope with 6V 20W Halogen Lamp A11.0303

- A11.0303 Microscope,
- Achromatic Objectives
- WF10x/18mm Wide Field Eyepiece
- Monocular, Binocular, Trinocular Head Optional
- Mechanical stage

Specification	A11.0303 Biological Microscope	M	D	B	T
Head	Monocular head, inclined 45°, rotatable 360°	•			
	Dual Viewing head, inclined 45°, rotatable 360°		•		
	Binocular head, inclined 30°, interpupillary distance 55-75mm			•	
	Trinocular head, inclined 30°, interpupillary distance 55-75mm				•
Eyepiece	WF10x/18mm				
Nosepiece	Quadruple nosepiece				
Objective	Achromatic Objective 4x/NA 0.1				
	Achromatic Objective 10x/NA 0.25				
	Achromatic Objective 40x(s)/NA 0.65				
	Achromatic Objective 100x(s,oil)/NA.1.25				
Stage	Mechanical stage size 140x120mm, moving range 70x40mm, minimum vernier division 0.1mm				
Focusing	Coaxial coarse & fine focusing knob, coarse focusing range 20mm, fine focusing division 0.002mm				
Condenser	N.A.1.25 Abbe condenser, up/down range 15mm				
Filter	Blue filter				
Light Source	6V 20W Halogen lamp, brightness adjustable				
Optional Accessory		Item No.			
Head	Seidentopf binocular head, inclined 30°, interpupillary distance 55-75mm	A53.0304-01			
	Seidentopf Trinocular head, inclined 30°, interpupillary distance 55-75mm	A53.0304-02			
Eyepiece	WF16x	A51.0301-16			
	WF20x	A51.0301-20			
Objective	Achromatic 20x (S)	A52.0301-20			
	Achromatic 60x (S)	A52.0301-60			
	Plan 4x	A52.0303-04			
	Plan 10x	A52.0303-10			
	Plan 20x(S)	A52.0303-20			
	Plan 40x(S)	A52.0303-40			
	Plan 60x(S)	A52.0303-60			
Condenser	Plan 100x(O,S)	A52.0303-100			
	Dark field condenser, dry	A57.0310-01			
	Dark field condenser, immersion	A57.0310-02			



Seidentopf Binocular Head is Optional

A11.0303-B



Seidentopf Triinocular Head is Optional

www.cnoec.com

A11.0303-T

www.cnoec.com

A11.0303-M

A11.0303-D



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China