



Binocular / Trinocular 1000x Fluorescence Microscope with Halogen Lamp A16.0301

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A16.0301
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Highlight: **inverted fluorescence microscope, epifluorescence microscope**



for more products please visit us on cnoec.com

Product Description

Binocular / Trinocular 1000x Fluorescence Microscope with Halogen Lamp A16.0301

- Fire-new ergonomic design, steady system structure, easy operation, is suitable for various working environments.
- Widely apply to clinical diagnosis, teaching experiment, pathological test and other micro-fields.

Specification	A16.0301 Metallurgical Microscope	B	T
Head	Binocular Head, inclined 30°, interpupillary distance 55-75mm	●	
	Trinocular Head, inclined 30°, interpupillary distance 55-75mm		●
Eyepiece	WF10x/18mm		
Nosepiece	Quadruple nosepiece		
Objective	Achromatic 4x/0.1		
	Achromatic 10x/0.25		
	Achromatic 40x/0.65(S)		
	Achromatic 100x/1.25(O,S)		
Stage	Mechanical stage size 140x140mm, moving range 70x50mm, minimum vernier division 0.1mm		
Condenser	N.A.1.25 Abbe condenser, up/down range 15mm		
Focusing	Coaxial coarse & fine focusing knob, coarse focusing range 20mm, fine focusing division 0.002mm		
Fluorescent Unit	A57.0302Epi-Fluorescent Unit		
	--B(Blue), G(Green), U(Ultraviolet), V(Violet) exciting light filter system		
	--100W super high pressure spherical mercury lamp --Ultraviolet light resist screen		
Filter	Yellow filter Φ43mm, Green filter Φ43mm		
Light Source	12V 50W Halogen lamp, brightness adjustable		
Optional Accessories		Item No.	
Eyepiece	WF16x	A51.0301-16	
	WF20x	A51.0301-20	
Objective	Achromatic 20x (S)	A52.0301-20	
	Achromatic 60x (S)	A52.0301-60	
	Plan 4x	A52.0303-04	
	Plan 10x	A52.0303-10	
	Plan 20x(S)	A52.0303-20	
	Plan 40x(S)	A52.0303-40	
	Plan 60x(S)	A52.0303-60	
	Plan 100x(O,S)	A52.0303-100	
Condenser	Dark field condenser, dry	A57.0310-01	
	Dark field condenser, immersion	A57.0310-02	
Adapter	CCD adapter	A55.0302	
	Photo adapter	A55.0301	
	Digital camera adapter	A55.0303	



A16.0301-B



A16.0301-T



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China