



Long Working Distance Stereo Optical Microscope with Screen Halogen Lamp A22.0302

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE, Rohs
- Model Number: A22.0302
- Minimum Order Quantity: 1
- Price: FOB \$100~1000, Depend on Order Quantity
- Packaging Details: Carton Packing
- Delivery Time: 5-30
- Payment Terms: L/C, T/T, Western Union
- Supply Ability: 10000



Product Specification

- Highlight: stereo dissecting microscope, stereo binocular microscopes



More Images



Product Description

Long Working Distance Stereo Optical Microscope with Screen Halogen Lamp A22.0302

- Original creation & designing in combination with visual project
- Easy changable objective 2x, 4x, 6x, 8x, 10x, humanisic operating system
- Suitable for long time working, comfortable operating and viewing improve efficiency
- Super long working distance, super large field of viewing & depth of field.
- Used in electrical industry, precision project, medical diagnosis and so on.

Our Product

A22.0302 Stereo Microscope				
Objective	Magnification	2x	4x	6x
	Working Distance (mm)	208	98	80
	Viewing Field (mm)	68	34	22.7
Base Plate Size	400x300x20mm			
Vertical arm height	500mm			
Travel	200mm			
Illumination	AC 110V/220V, 12V 20W halogen lamp with cup (2 pcs), Outside LED light adjuster			
A12.0205 Laboratory Biological Microscope Optional Accessories				
Objective	Magnification	Working Distance (mm)	Viewing Field (mm)	Item No.
	8x	58	17	A52.0380-08
	10x	46	13.6	A52.0380-10



E22.0302



E22.0302



E22.0302



Large View Field Glass



Objective



Observing a Screw



Image Effect



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China