



## 7x-45x Binocular / Trinocular Jewelry Microscope With Zoom Ratio 1 / 9 A24.0301

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A24.0301
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Research
- Certification: CE
- Drawtube: Binocular
- Theory: Jewelry Microscope
- Highlight: **diamond microscope, gemological microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### 7x-45x Binocular / Trinocular Jewelry Microscope With Zoom Ratio1 / 9 A24.0301

A24.0301 Gem Microscope, with high definition and high resolution rate of solid image, providing special clamp to hold gem conveniently, the inside crystal structure of gem can be shown clearly under 6000K light source, it is very suitable for professional gem and diamond checking field.

Specification	A24.0301	A	B
Head	Binocular head,interpupillary distance 55-75mm	●	
	Trinocular head,interpupillary distance 55-75mm		●
Total Magnification	7x-45x		
Eyepiece	Super WF10x/23mm eyepiece		
Objectives	Zoom objective 0.7x-4.5x		
Body	Body can be inclined 45° for comfortable operating		
Gem clamp	Alligapop clamp, steel clamp, plat clamp included		
Condenser	Auxiliary lens optional : 0.75x (W.D.105mm), 1.5x(W.D.47mm), 2x(33mm)		
Focus	Focusing range 125mm, working distance 95mm		
Illumination	Bright/Dark field illuminator, with iris adjustable 2-41mm, 8W 6000K fluorescent lamp		



A24.0301-A



www.cnoec.com



A24.0301-B



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China