



## Infinity Compensation Wide Field Microscope A12.0302 Binocular Head , Y Style Body

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.0302
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Feature: High Resolution
- Type: A12.0302
- Standard: CE, Rohs
- Original: China
- Shape: "Y" style Design
- Quadruple: Nosepiece
- Highlight: **binocular compound microscope, optical compound microscope**



### More Images



## Product Description

### Infinity Compensation Free Binocular Head Compound Optical Microscope A12.0302

- CCIS Infinity Chromatic Aberration Revise System, High resolution image.
- Improve in optical capability, wide-field, better image quality.
- Main body of "Y" style design, stable & reliable whole structure.
- Humanized configuration design, lower draw operation, makes the user more comfortable.
- Knocked-down type and double layers mechanical stage.

	A12.0302-B	A12.0302-T
Head	Compensation Free Binocular Head	Compensation Free Trinocular Head
	Interpupillary Adjustable Distance: 50-75mm	
Eye piece	WF10X (Viewing Field Diameter: $\varnothing$ 20mm)	
Nose piece	Quadruple Nosepiece	
Objective	Infinity Plan Achromatic Objectives: 4x $\infty$ /NA 0.1 10x $\infty$ /NA 0.25 40x(s) $\infty$ /NA 0.65 100x(s, oil) $\infty$ /NA.1.25	
Working Stage	Mechanical stage size 185x145mm, moving range 75x50mm, minimum vernier division 0.1mm	
Focusing	Coaxial Coarse & Fine Focus Knob, With Tension Adjustable And Limit Stopper, Minimum Division 0.002mm, Moving Range 0~20mm	
Condenser	N.A.1.25 Abbe condenser	
Light Source	6V 30W Halogen lamp	






**Opto-Edu (Beijing) Co., Ltd.**

---

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China