



10x - 80x Binocular Jewelry Microscope Zoom Ratio 1 / 4 A24.0401

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A24.0401
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Research
- Certification: CE
- Drawtube: Binocular
- Theory: Jewelry Microscope
- Magnification: 10x - 80x
- Illumination: Halogen Darkfield/Brightfield
- Highlight: **laboratory binocular microscope,
diamond microscope**

for more products please visit us on cnoec.com

Product Description

10x - 80x Binocular Jewelry Microscope Zoom Ratio 1 / 4 A24.0401

Gemological microscope is designed for jewelers and gem stone experts, to be one of the most important tool in their jobs. A24.0401 gemological microscope is especially designed to view precious stone samples and the pieces of jewelry that contain them.

A24.0401 Gem Microscope	
Head	Binocular head, interpupillary distance: 55-75mm, diopter adjuster: ± 5 diopter
Body	Optic head and stand tilt to a more natural 40° viewing angle
Eyepiece	WF10x
Zoom Ratio	1:4
Objective Zoom Range	10x-80x
Focusing	Focus range: 50mm
Working Distance	95mm
Illumination	Voltage input 90-230V
	Adjustable iris diaphragm to control width of 2-42mm light beam
	Built-in optical-fiber illumination system
	Tungsten-halogen darkfield/brightfield illumination



A24.0401



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China