



A43.1904 OPTO EDU Binocular Edan Colposcope Continuous Zoom 0.4x-2.75x

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A43.1904
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Observation Angle: 45° Tilted Binocular Observation Unit
- Eyepiece: 10x
- Interpupillary Distance: 50mm-80mm
- Diopter: $\pm 5D$
- Objective Focal Length: F=300mm (M58x0.75mm)
- Magnification Changer: 0.4x-2.75x
- Total Magnification: 2.4x-16x
- Linear Field: 80mm-13mm
- Fine Focusing Distance: 20mm
- Floor Stand & Arm: Universal Stand,
- Highlight: **Binocular Edan Colposcope,
Continuous Zoom Edan Colposcope,
2.75x Edan Colposcope**

Product Description

A43.1904 OPTO EDU Binocular Edan Colposcope Continuous Zoom 0.4x-2.75x

Table Of Specification For A43.1904 Series Colposcope

Observation Angle	45° Tilted Binocular Observation Unit
Eyepiece	10x
Interpupillary Distance	50mm-80mm
Diopter	±5D
Objective Focal Length	F=300mm (M58x0.75mm)
Magnification Changer	0.4x-2.75x
Total Magnification	2.4x-16x
Linear Field	80mm-13mm
Beam Splitter	Included, 50:50
Video Camera Adapter	C-mount Included, For 1/3" Digital Camera
Fine Focusing Distance	20mm
Illumination System	Coaxial Illumination, 10W LED, Brightness Adjustable, Illumination>30000lx
Power Supply	Power Supply Input: AC110~240V, 50/60Hz, With Power Supply Adapter Output DC12V3A, 36W,
Filter	Built-in Green And Blue Filters
Floor Stand & Arm	Universal Stand, Balancing Arm With Counterweight & Lock Short Column Fitted On Five-star Base With Casters
Optional Accessories	Demonstrator, Camera (A59.4910, A59.4231), SLR Camera Adapter, Monitor, Objective (F=230mm & F=280mm)



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085


A43.1904

Colposcope, Continuous Zoom 0.4x-2.75x, Universal Stand



A43.1904 Optional Accessories

OPTO-EDU





A43.1904 With Different Digital Camera (Optional) + Monitor/Computer (Optional)



A43.1904



Binocular Head 45° Tilted, IPD50-80mm, Diopter ±5D, Eyepiece WF16x, Continuous Zoom 0.4x~2.75x, Objective F=300mm



A43.1904 With Different Digital Camera (Optional) + Monitor/Computer (Optional)



Opto-Edu (Beijing) Co., Ltd.

0086 13911110627

sale@optoedu.com

cnoc.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China