



Opto Edu A14.0301 Inverted Optical Microscope Binocular 400x Abbe Condenser 12V 50W

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A14.0301
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month

OPTO-EDU



Product Specification

- Eyepiece: WF10x/18mm
- Nosepiece: Quadruple Nosepiece
- Objectives: Long Working Distance Plan Achromatic
- Working Stage: Mechanical Stage Size 190x145mm, Moving Range 75x50mm, Minimum Vernier Division 0.1mm
- Condenser: N.A.1.25 Abbe Condenser, Up/down Range 15mm
- Light Source: 12V 50W Halogen Lamp, Brightness Adjustable
- Highlight: **Binocular Inverted Optical Microscope, 400x Inverted Optical Microscope, 12V Inverted Light Microscope**

OPTO-EDU



More Images

OPTO-EDU



OPTO-EDU



Product Description

Our Product Introduction

Specification	A14.0301 Inverted Microscope	B	T
Head	Binocular Head, inclined 30°, interpupillary distance 55-75mm	.	
	Trinocular Head, inclined 30°, interpupillary distance 55-75mm		.
Eyepiece	WF10x/18mm		
	WF16x/11mm		
Nosepiece	Quadruple nosepiece		
Objective	Long Working Distance Plan Achromatic 4x/0.1, W.D.=17.4mm		
	Long Working Distance Plan Achromatic 10x/0.25, W.D.=8.0mm		
	Long Working Distance Plan Achromatic 25x/0.40, W.D.=5.0mm		
	Long Working Distance Plan Achromatic 40x/0.6(S), W.D.=3.0mm		
	Auxiliary objective 1.25x Thickness in the bottom of culature ware 1.5mm		
Stage	Mechanical stage size 190x145mm, moving range 75x50mm, minimum vernier division 0.1mm		
Condenser	N.A.1.25 Abbe condenser, up/down range 15mm		
Focusing	Coaxial coarse & fine focusing knob, coarse focusing range 30mm, fine focusing division 0.002mm		
Filter	Yellow filter Φ43mm, Green filter Φ43mm		
Light Source	12V 50W Halogen lamp, brightness adjustable		
Optional Accessories		Item No.	
Adapter	CCD adapter	A55.0302	
	Photo adapter	A55.0301	
	Digital camera adapter	A55.0303	



OPTO-EDU (BEIJING) CO., LTD.

OPTO-EDU

F-1501 Wanda Plaza, No.18 Shijingshan Road, Beijing 100043, China
Tel:+8610 88696020 Fax:+8610 88696085


A14.0301 Inverted Microscope



for more products please visit us on cnoec.com

Sliding Head

130°, Interpupillary Distance 55-75mm
With WF10x/18mm Eyepiece



Mechanical Stage

Size 190x145mm, Moving Range 75x50mm, Minimum Vernier Division 0.1mm



12V 50W Halogen Lamp

brightness Adjustable
Filters: Yellow Filter $\Phi 43\text{mm}$, Green Filter $\Phi 43\text{mm}$




Long Working
Distance Plan
Achromatic
Objectives:
10x,25x,40x

A14.0301 Specification

OPTO-EDU





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China