



Nikon Series Microscope Accessories , A55.2207 C-Mount Adapter

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC
- Certification: CE
- Model Number: A55.2207
- Minimum Order Quantity: 1
- Packaging Details: Carton
- Supply Ability: 1000



Product Specification

- Highlight: microscope accessory, parts of the microscope



for more products please visit us on cnoec.com

Product Description

Nikon Series Microscope Accessories , A55.2207 C-Mount Adapter

A55.2207, Nikon C Mount, The Low Magnification C-mount, U-TV Series is low-magnification TV adapter with C-mount which allows a digital imaging device such as a digital TV camera, CCD camera, UHCCD camera, or CMOS camera to capture wide-field images through microscope trinocular head. Use of telecentric optics helps reduce the occurrence of light deficiency in the peripheral sections. The transmittance at the infrared frequency band is extended to 1,000 nm.

- Low Magnification C-mount Adaptor for Wide Field Image
- Adjustable Focusing or Fixed Lens Inside
- Telecentric Optics with Low Light Deficiency
- Parfocal with Different Microscope Objective Lenses
- Can Be Installed in UIS Trinocular Tube
- Diffraction Limited MTF Performance

A55.2207	0.35X	0.55X	0.70X	1.0x
Magnification	0.35x	0.55x	0.70x	1.0x
Field Size (mm)	22x0.35=7.7	22x0.55=12.1	22x0.7=15.4	22x1.0=22
Compatible Sensor	1/2.5", 1/3", 1/4"	1/1.8", 1/2", 1/2.5"	2/3", 1/1.8", 1/2"	1", 2/3"



A59.2207 C-Mount, Nikon Series






Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China