



144 LEDs LED Ring Light Microscope Accessories Adjustable 64 LEDs A56.1211

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.1211
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE|Rohs
- Product Name: Microscope Accessories
- Color Temp: 3500k
- Power Input: 90-265VAC
- Power Output: 22VDC
- Highlight: microscope accessory, parts of the microscope



More Images



for more products please visit us on cnoec.com

Product Description

144 LEDs LED Ring Light Microscope Accessories Adjustable 64 LEDs A56.1211

Description:

- Long Life
- Light Adjustable
- Fits zoom microscope

Model	A56.1211-B48	A56.1211-B56	A56.1211-B64	A56.1211-B144
LED Bulbs	48 pcs	56 pcs	64 pcs	144 pcs
Output	12V / 0.5A			
Brightness Adjustment	0-100% Adjustable			
Working Distance	50-250mm			40-250mm
Brightness	Around 12000lux			Around 18000lux
Working Life	10000 Hours			



A56.1211 LED Ring Light



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

