



Mobile Phone Microscope Accessories , Micorscope A55.9010 Universal Holder

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A55.9010
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE|Rohs
- Product Name: Microscope Accessories
- Magnification: 3.8x - 230x
- Working Distance: 101mm
- View Field: Dia. 2.3~62mm
- Highlight: **microscope accessory, compound microscope parts**



for more products please visit us on cnoec.com

More Images



Product Description

Mobile Phone Microscope Accessories , Micorscope A55.9010 Universal Holder

we have our own injection machines, we are a manufacturer, can also accept OEM order. Or if you have you own idea on new style products and need a factory to cooperate, we promise we should be the most suitable factory you need

Mobile Phone Width 55~98mm, Eyepiece Dia26.4~46.4mm of All Kinks Microscope, Binocular, Telescope Connection



A55.9010



A55.9010 Use with Mobil Phone






A55.9010 Use With Microscope




Observe Effect Images:



A55.9010 Use with Telescope






Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China