



8x - 50x Trinocular Stereo Zoom Microscope with Track Stand LED Light A23.1501

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 8x-50x
- Applications: Industry|Gem|Medical|Clinic
- Certification: CE|Rohs
- Drawtube: Binocular|Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo dissecting microscope, stereo binocular microscopes**



More Images



Product Description

8x - 50x Trinocular Stereo Zoom Microscope with Track Stand LED Light A23.1501

A23.1501 Series zoom stereo microscope:

- Well-coated optical system can provide sharp & clear images with extremely flatness and contrast.
- Ergonomic design can make users feel comfortable after long-time operation
- Standard zoom magnification 8X-50 with zoom ratio 1:6.3. Standard working distance is 115mm. With additional objectives, it can be increased to 211mm.
- Configured with WF10X/22mm eyepiece. Each eyepiece allows diopter adjustment to suit different users.
- Full set of eyepiece, auxiliary objectives and other accessories can fulfill different needs.

A23.1501 Zoom Stereo Microscope		B	B	B	B	T	T	T	T	
		1	2	3	4	1	2	3	4	
Head	Binocular head, 45° inclined, 360° rotatable, Interpupillary range from 52mm-75mm	o	o	o	o					
	Trinocular head, 45° inclined, 360° rotatable, Interpupillary range from 52mm-75mm . Distribution is 3:7 for trinocular					o	o	o	o	
Magnification	Standard Zoom Magnification: 8x-50x. With additional eyepiece & objectives, can reach 4x-200x.									
Eyepiece	WF10x/22mm									
Working Distance	Standard Working Distance is 115mm. With additional objectives can be up to 45mm-211mm.									
Auxiliary Objective	0.5x/210mm; 0.7x/130mm; 2x/45mm									
Focusing	Coarse Focusing Knob with tension adjustable. Adjusting range 100mm									
Base	Track Big Fan Stand, No Light	o				o				
	Track Big Fan Stand, Upper/Bottom LED		o				o			
	Track Big Square Stand, No Light			o				o		
	Track Big Square Stand, Upper/Bottom LED				o				o	
Optional Accessories									Item No.	
Eyepiece	WF10x/22mm, Diopter adjustment D-8~+5								A51.1522	-10
	WF15x/16mm, Diopter adjustment D-8~+5								A51.1522	-15
	WF20x/12.5mm, Diopter adjustment D-8~+5								A51.1522	-20
Light	LED Ring Light , 60pcs								A56.1522	



www.cnoec.com

A23.1503-B4



Accessories:



Standard Eyepiece
WF10x/22

Optional Eyepiece:
WF10x/22, WF15x/16, WF20x/12 with Diopter Adjustment



A56.1522
Optional LED Ring Light 60pcs

A56.1522 Use with A23.1501 Microscope



Auxiliary Objectives: 0.5x/220mm; 0.7x/137mm; 2x/45mm

Product Size:






Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China