



Binocular Laboratory Compound Microscope 40 - 1000x with 3W LED A12.1503

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1503
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Binocular
- Illumination: LED
- Theory: Biological Microscope
- Optical System: Finite, 160mm
- Highlight: **binocular compound microscope ,
wide field microscope**

for more products please visit us on cnoec.com

Product Description

Binocular Laboratory Compound Microscope 40 - 1000x with 3W LED A12.1503

- Binocular And Trinocular Head Available, Interpupillary Distance 48~75mm
- Both Eyepiece Tube Diopter Adjustable
- Wide Field Eyepiece WF10X/18mm
- 3W/LED

Specific ation	A12.1503-A	A12.1503-B
Optical System	Finite, 160mm	
Head	Seidentopf Binocular head,30° inclined, 360° rotatable, interpupillary distance 48-75mm,	Seidentopf Trinocular head,30° inclined, 360° rotatable, interpupillary distance 48-75mm,
Eyepiece	WF10x/18mm	
Objective	Achromatic objectives 4x,10x,40x & 100x (oil)	
Nosepiece	Forward Quadruple Nosepiece	
Focus System	Coaxial Coarse and Fine Adjustment, Fine Division 0.002mm, Coarse Stroke36mm	
Stage	Double Layer Mechanical Stage 142x132mm Moving Range75x50mm,ow position and right hand controls, X-Y coaxial control knob	
Condenser	Abbe N.A.1.25 with iris diaphragm	
Light Source	External Illumination, 3W/LED or 6V20W Halogen Lamp	
Adapter	0.5x, C-Mount	
A12.1503 Laboratory Microscope Optional Accessories		Item No.
Eyepiece	WF16	A51.1502-16



A12.1503



Forward Quadruple Nosepiece



Abbe N.A. 1.25 with Iris Diaphragm



Infront View of A12.1503



Side View of A12.1503



Opto-Edu (Beijing) Co., Ltd.

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China